



ACCELERADOR

Descripción:

Líquido acelerante del tiempo de fraguado y endurecedor para yesos de clase III

Utilización:

Se presenta en forma concentrada para ser diluido con agua desmineralizada en una proporción del 10 al 50 % en volumen.

Método de trabajo:

Diluir el ACCELERANTE PARA YESOS en la proporción deseada según la disminución del tiempo de trabajo que se necesite.

Ejemplo de utilización con un yeso tipo DURGUIX (proporción de mezcla 28 ml / 100 g)

Cantidad de ACCELERANTE PARA YESOS (ml)	3	7	14
Cantidad de agua desmineralizada (ml)	25	21	14
Dilución	10%	25%	50%
Disminución del tiempo de fraguado (min)	0.5	2	3.5

Presentación:

1195-100 Botella de 1 l

Precauciones de uso:

La estabilidad del producto es de 18 meses (respetando las instrucciones de conservación.)

Estas indicaciones están basadas en nuestro nivel de conocimiento actual del producto basado en nuestras propias experiencias y ensayos y deben ser comprendidas sólo como valores indicativos no entendiéndose en ningún caso como una garantía del producto.



protechno

Polígono Empordà Internacional
17469 VILAMALLA, GIRONA, SPAIN

Tel.: (34) 972 52 61 69

protechno@protechno.com

www.protechno.com



ACCELERADOR

Description:

Accelerating liquid to reduce the setting time and increase the hardness of type III gypsum

Uses

It's a concentrate liquid to dilute with distilled water. Dilution from 10 to 50 % in volume

Working instructions

Dilute ACELERADOR in the wished proportion to achieve the needed decrease in working time

Use example with DURGUIX (type III gypsum. Mixing proportion de 28 ml / 100 g)

Amount of ACELERADOR (ml)	3	7	14
Amount of distilled water (ml)	25	21	14
Dilution	10%	25%	50%
Decrease in setting time (min)	0.5	2	3.5

Supply

1195-100 1 l Bottle

Use cautions:

The stability of the product is 18 months respecting the storage.

The instructions that appear on this sheet are based on our current knowledge and are given only as information without that implies a guarantee from Protechno.



protechno

Polígono Empordà Internacional

17469 VILAMALLA, GIRONA, SPAIN

Tel.: (34) 972 52 61 69

protechno@protechno.com

www.protechno.com