ORTOGUIX PLASTER

Descripción

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ORTOGUIX PLASTER

Description

Special plaster for study models in orthodontics. Type II, high quality ISO 6873

Uses

Special plaster for study in orthodontics. Easy to cut

Advantages

This product satisfies the highest requirements of the professionals.

Good expansion

Very good definition of details.

Snow white colour that doesn't get yellow with time

Physical characteristics

Working time:	7 min.	
Initial setting:	10 min. (Vicat test according to ISO 6873)	
Final setting :	13 min. (Vicat test according to ISO 6873)	
Setting time:	30 min.	
Setting expansion (2h):	0.16 %	
Hardness (24h) :	95 N / mm ²	
Compressive strength (1h) :	25 N / mm ²	

Working instructions

Powder/water ratio:	100 g / 37-38 ml
Working temperature:	20-23°C
Mixing:	mechanical under vacuum: 30 sec. manual: 60 sec.

Add the powder to the demineralised water following the recommended ratio. The initial consistency of the slurry must be thick. In this moment don't add more water above the recommended.

While the mixing time is complete, the product acquires the suitable fluidity and stability. When this mixing time has finished, the mix must have a soft and high tyxotropic consistency, flowing easily over the model with a low vibration.

After 4 minutes, stop the vibration to guarantee a good crystallisation of the product.

Cautions

1.- USE CAUTIONS:

Use always-clean instruments

Do not mix the product in a container where investments have been mixed.

Use always demineralised water to do the mix. The water from the public mains can affect the correct setting process of the gypsum, making it faster or slower

2.- SAFETY CAUTIONS:

Do not inhale the powder. It's recommended the use of a protective mask.

Keep away from children.

Handle the product in a well-ventilated place.

The handling of this product has to be done only by qualified people

In the event of ingestion, consult a doctor.

For further information, see the safety data sheet for the product.

Storage and expiry

The stability of the product is 18-36 months depending on the package and respecting the storage conditions. (Note: The expiry date is considered with the product in his original packing and before been opened).

Storage: Keep the packages tightly closed in a dry place away from moisture.

After each use the package must be closed to prevent the moisture absorption.

Supply

White

1151-200 Pail 20 Kg

Additional information

The instructions that appear on this sheet are based on our current knowledge and are given only as information without that implies a guarantee from Protechno.

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