## RESIN GIPS

NEW IMPROVED  
FORMULA

### Description

Extra hard high density and resin die-stone. Type IV, high quality ISO 6873

### Uses

Ideal for models in skelectric prostheses, crowns, bridges, inlays, combined works and master models with the technique of precious metals and ceramic .

### Advantages

New improved formula:

- bigger resin content
- greater resistance to abrasion

Due to the formulation with resin, you obtain very hard and even surfaces

Compatible with all kind of impression materials

Low setting expansion which provokes a good dimensional stability

Easy to cut avoiding the formation of dust in the laboratory

Very fluid gypsum : easy to use.

Exact reproducibility.

Snap setting effect

### Physical characteristics

Working time: 7 min

Initial setting: 9 min ( Vicat test according to ISO 6873)

Final setting : 10 min (Vicat test according to ISO 6873)

Setting time: 30 min

Setting expansion (2h): 0.10%

Compressive strength (1 h): 50 MPa

### Working instructions

Powder/water ratio: 100 g / 20 ml

Working temperature: 20-23°C

Mixing: mechanical under vacuum: 30 sec.  
manual: 60 sec.

Add the powder to the demineralised water following the recommended ratio.

**Warning:** At the beginning of the mixture may appear that more water is needed, but do not add more. Once the mixing time has finished, the right consistency will be achieved.

**It is highly recommended the vacuum mixing**

When this mixing time has finished, the mix must have a very fluid consistency, flowing easily over the model with a low vibration. After 4 minutes, stop the vibration to guarantee a good crystallisation of the product.

### Cautions

#### 1.- USE CAUTIONS:

Use always-clean instruments

Do not mix the product in a container where investments have been mixed.

Use always demineralised water to do the mix. The water from the public mains can affect the correct setting process of the gypsum, making it faster or slower

#### 2.- SAFETY CAUTIONS:

Do not inhale the powder. It's recommended the use of a protective mask.

Keep away from children.

Handle the product in a well-ventilated place.

The handling of this product has to be done only by qualified people

In the event of ingestion, consult a doctor.

For further information, see the safety data sheet for the product.

### Storage and expiry

The stability of the product is 36 months respecting the storage conditions. (Note: The expiry date is considered with the product in his original packing and before been opened).

Storage: Keep the packages tightly closed in a dry place away from moisture.

After each use the package must be closed to prevent the moisture absorption.

### Supply

#### Colour apricot

1107-120 Eco – pack: 12 Kg ( 6 x 2 Kg Al bag)

1107-k20 Test-pack: 2 Kg (1 x 2 Kg)

#### Additional information

The instructions that appear on this sheet are based on our current knowledge and are given only as information without that implies a guarantee from Protechno.

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