TECHNOJEL-PRO

Descripción

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TECHNOJEL-PRO

Description

High hardness reusable duplicating material, Agar-Agar based.

Use

TECHNOJEL-PRO is specially designed for the production of dentures with the acrylic pouring technique. We recommend the use of ECO-CRYL PRO resin

It can also be used to produce phosphate or polysilicate investment duplicates within the technique of reproduction of models.

Do not use to duplicate plaster.

Advantages

Due to its improved hardness, with TECHNOJEL-PRO you will be able to obtain duplicates with high reliability and great accuracy and detail definition.

Hardens with no bubbles.

Cost reduction due to continuous use and higher output.

Easy to use: fusion with duplicating machine, boiling water bath or microwave

Working instructions

MELTING WITH DUPLICATING MACHINE OR BOILING WATER BATH

- Cut the gelatine in small pieces (about 4 cm) and introduce them in the duplicating unit or boiling water bath
- Melt the product with continuous stirring. Do not exceed 90 °C, the gel must not boil. Check that no small pieces are kept unmelted.
- Leave the product cool down to 45 – 50 °C (optimum working temperature). At this temperature the product is ready to be poured in the muffle.
- Soak the master model before the duplication to improve the process
- Gelation time is approximately 30 minutes in a metallic muffle.
- In plastic muffle this time is substantially prolonged.

MICROWAVE MELTING

- Cut the gelatine in small pieces (about 4 cm) and put them into an appropriate plastic container
- 500 g TECHNOJEL-PRO: heat 4 minutes at 800 W, stir carefully and heat again 1 minute at 400 W
- For bigger quantities, extend the heating time, do not increase the power
- The gel must not boil. Check that no small pieces are kept unmelted.
- Leave the product cool down to 40 – 50 °C . Stir occasionally
- Once the working temperature is achieved, the product is ready to be poured in the appropriate muffle.
- Gelation time is approximately 40 minutes.

Duplicating:

- After setting, open the flask and remove the model carefully
- If the gelatine is too cold, it can lead to a longer setting time and therefore a weaker and rougher duplicate.
- Pour the investment or acrylic resin following the manufacturer working instructions. We recommend our acrylic resin ECO-CRYL PRO for the production of dentures with the pouring technique and CALIBRA-M to obtain accurate models for the framework technique .

Re-use:

This product can be re-used following these indications:

- Wash the used gelatine carefully removing all investment remains. Keep it in a closed container.
- Cut in small pieces and melt again with the duplicating machine or microwave as if it was a new product
- If the gelatine it's too thick, small quantities of distillate water (lost by evaporation) can be added

When the gelatine has been re-used several times, it can lose its properties (fluidity, plasticity, colour, ...) If so, discard it and use a brand new product.

TECHNOJEL-PRO can be re-used about 6 times without losing its properties only if instructions are followed

Do not melt used and new material together

Physical characteristics

Melting point:	90 °C.	(194°F)
Solidification point:	38 °C app.	(82°F)
Use temperature:	45 - 50 °C	(113 - 122°F)
Gelation time:	30 – 40 min	

Availability

5512-060	6 Kg container
5512-010	1 Kg container

Auxiliary products for the pouring technique:

5555-100	Pouring flask
9045-121	ECO-CRYL PRO resin . Kit Standard

Expiry and Storage

The stability of the product is 2 years (respecting the conservation instructions.)

The product should be conserved inside its packing properly closed and in a dry atmosphere.

Additional information

The information given in this sheet is based on our level of knowledge at the moment at the date of printing and is not intended to be a warranty for the product.

Polígono Empordà Internacional
17469 VILAMALLA, GIRONA , SPAIN
Tel.: (34) 972 52 61 69
protechno@protechno.com
www.protechno.com

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