## POLICROM-NEUTRALIZADOR

### Descripción

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## POLICROM-NEUTRALIZADOR

### Description

POLICROM-NEUTRALIZADOR is a product to neutralise electrolytic polishing liquids (POLICROM), specially designed to solve the disposal problems of corrosive products without environmental and personnel dangers.

### Instructions for use

The following steps should be performed to achieve a correct neutralisation:

Mix always POLICROM-NEUTRALIZADOR with the electrolytic liquid using the original container.

It is recommended to add small amounts of used electrolytic liquid and shake the mixture every addition.

Do not add the total final amount all at once; a sudden temperature increase of the mixture could happen.

Once the container is full, the liquid will be neutralised and will be no longer corrosive.

The obtained mixture should be disposed according to the local regulation. It still contains heavy metals in solution which come from the metal polishing.

### Further information

The optimum mixing proportion is 1:1.  
Lower or higher proportions will not be valid for total neutralization.

pH of the final mixture neutralizador /electrolytic liquid is about 10-10,5 (at a Temperature of 20°C).

Dispose only when the container is full and the mixture is cold.

Do not use with electrolytic liquid of another manufacturer.

### Packaging

7891-050 5 L bottle ( bottle capacity : 10 L )

### Expiry and storage

The stability of the product is 3 years respecting the storage conditions. (Note: The expiry date is considered with the product in its original packaging and before been opened).


Keep the product in the original container in a dry place away from moisture.

Store in a safe place, the product is corrosive.

After each use the package must be closed to prevent the moisture absorption.

### Precautions

Corrosive

### Additional information

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet. Order them to your local dealer or visit our website to download them.

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