TECHNOSIL SUPER

Descripción

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TECHNOSIL SUPER

Description

Separating fluid for acrylic resins

Composition

This is an alginate-based separating media with emulsified oils.

Application

Can be used on the cold or preheated model, both in the flasking technique (removable prosthesis, crowns and bridges) and in hydro pneumatic polymerisation.

It can be applied several times on any kind of gypsum without danger of peeling off.

The emulsion forms a fine film that creates an effective barrier between the gypsum and the resin.

We recommend pouring the necessary amount into the dispenser and spreading with a soft brush; the operation should be repeated until the gypsum has a uniformly glossy appearance.

Repeating the operation for a second time is usually sufficient, although this will depend on the degree of absorption of the various gypsums used.

When the product has remained unused for a long time some components may have separated out.

THE PRODUCT SHOULD ALWAYS BE STIRRED BEFORE EACH USE.

Storage and stability

The stability of the product is 2 years respecting the storage conditions. (Note: The expiry date is considered with the product in his original packing and before been opened).

Storage: Keep the product tightly closed in its original package and in a dry cool place.

Packing

3211-455 Bottle 450 ml

3211-M55 Sample Bottle 50 ml

Additional information

The information given in this sheet is based on our level of knowledge at the moment at the date of printing and is not intended to be a warranty for the product.

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