MUFLA-GIPS

Description

Fluid plaster for flasking. Type II, ISO 6873

Uses

Very fluid plaster to be poured into the upper half of the flask without vibration.

Recommended for all kind of flasks to make the upper half of the flask.

Advantages

Once set, the obtained surface is exceptionally smooth, which greatly reduces air pockets and bubbles in the acrylics.

Free of bubbles even without vibration and with manual mixing.

Does not adhere to acrylic denture or teeth

Natural white colour.

Physical characteristics

Working time:	3 – 4 min.
Initial setting:	5 min. (Vicat test according to ISO 6873)
Final setting :	8 min. (Vicat test according to ISO 6873)
Setting time:	25 min.
Setting expansion (2h):	0.30 %
Compressive strength (1h) :	15 N / mm ²

Working instructions

Powder/water ratio:	100 g / 35 – 40 ml
Working temperature:	20-23°C
Mix by hand	60 sec.

Add the powder to the water following the recommended ratio. When this mixing time has finished, the mix has a great fluid consistency, flowing easily into the flask even without vibration.

Cautions

1.- USE CAUTIONS:

Use always-clean instruments

Do not mix the product in a container where investments have been mixed.

Use always demineralised water to do the mix. The water from the public mains can affect the correct setting process of the gypsum, making it faster or slower

2.- SAFETY CAUTIONS:

Do not inhale the powder. It's recommended the use of a protective mask.

Keep away from children.

Handle the product in a well-ventilated place.

The handling of this product has to be done only by qualified people. In the event of ingestion, consult a doctor.

For further information, see the safety data sheet for the product.

Storage and expiry

The stability of the product is 24 months respecting the storage conditions. (Note: The expiry date is considered with the product in his original packing and before been opened).

Storage: Keep the packages tightly closed in a dry place away from moisture.

After each use the package must be closed to prevent the moisture absorption.

Supply

Colour: white

1184-200

Box 20 Kg

Additional information

The instructions that appear on this sheet are based on our current knowledge and are given only as information without that implies a guarantee from Protechno.

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