RESIN GIPS -Plus

Description

Extra hard resin added die-stone. Type IV, high quality ISO 6873

Uses

Specially designed for models in frame-work prostheses, crowns, bridges, inlays, combined works and master models with the technique of precious metals and ceramic.

Advantages

New improved formula:

- **higher resin content**
- **greater resistance to abrasion**

Due to the formulation with resin, you can obtain very hard and even surfaces

Compatible with all kind of impression materials

Controlled setting expansion which provides a good dimensional stability

Easy to cut avoiding the formation of dust in the laboratory

Very fluid gypsum: easy to use.

Exact reproducibility.

Snap setting effect

Physical characteristics

| | |
|-------------------------------|--|
| Working time: | 7 min |
| Initial setting: | 8 min (Vicat test according to ISO 6873) |
| Final setting: | 10 min (Vicat test according to ISO 6873) |
| Setting time: | 30 min |
| Setting expansion (2h): | 0.10 % |
| Compressive strength (1 h): | 50 MPa |

Working instructions

| | |
|----------------------|----------------------------------|
| Powder/water ratio: | 100 g / 20 ml |
| Working temperature: | 20-23°C |
| Mixing: | mechanical under vacuum: 30 sec. |

Add the powder to the demineralised water following the recommended ratio.

Warning: At the beginning of the mixing may appear that more water is needed, but do not add more. Once the mixing time has finished, the right consistency will be achieved.

It is highly recommended the vacuum mixing

When this mixing time has finished, the mix must have a very fluid consistency, flowing easily over the model with a low vibration. After 4 minutes, stop the vibration to guarantee a good crystallisation of the product.

Cautions

1.- USE CAUTIONS:

Use always-clean instruments

Do not mix the product in a container where investments have been mixed.

Use always demineralised water to do the mix. The water from the public mains can affect the correct setting process of the gypsum, making it faster or slower

2.- SAFETY CAUTIONS:

Do not inhale the powder. It's recommended the use of a protective mask.

Keep away from children.

Handle the product in a well-ventilated place.

The handling of this product has to be done only by qualified people

In the event of ingestion, consult a doctor.

For further information, see the safety data sheet for the product.

Storage and expiry

The stability of the product is 36 months respecting the storage conditions. (Note: The expiry date is considered with the product in his original packing and before been opened).

Storage: Keep the packages tightly closed in a dry place away from moisture.

After each use the package must be closed to prevent the moisture absorption.

Supply

Colour: grey

1108-062 Eco – pack: 6 Kg (3 x 2 Kg Al bag)

1108-K20 Test-pack: 1 Kg (1 x 2 Kg)

Additional information

The instructions that appear on this sheet are based on our current knowledge and are given only as information without that implies a guarantee from Protechno.

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