ROCANIT 0.08

Descripción

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ROCANIT 0.08

Description

Improved natural stone, low expansion. Type III, ISO 6873

Uses

Ideal for general laboratory uses, antagonist models, model bases with stumps, in the formation of impression trays and for any work that needs high precision and moderate hardness.

Advantages

This product satisfies the highest requirements of the professionals.

Very low setting expansion

Hard and even surfaces without bubbles

Good definition of details.

Very tyxotropic gypsum with good consistence for an easy handling

Physical characteristics

Working time: 5 – 6 min.

Initial setting: 7' – 7' 30" (Vicat test according to ISO 6873)

Final setting: 8' 30" – 9' (Vicat test according to ISO 6873)

Setting time: 20 min.

Setting expansion (2h): 0.08 %

Hardness (24h) : 137 N / mm²

Compressive strength (1 h) : 32 N / mm²

Working instructions

Powder/water ratio: 100 g / 30 ml

Working temperature: 20-23°C

Mixing: mechanical under vacuum: 30 sec.
manual: 60 sec.

Add the powder to the demineralised water following the recommended ratio. The initial consistency of the slurry must be thick. In this moment don't add more water above the recommended. While the mixing time is complete, the product acquires the suitable fluidity and stability.

When this mixing time has finished, the mix must have a soft and high tyxotropic consistency, flowing easily over the model with a low vibration. After 4 minutes, stop the vibration to guarantee a good crystallisation of the product.

Cautions

1.- USE CAUTIONS:

Use always-clean instruments

Do not mix the product in a container where investments have been mixed.

Use always demineralised water to do the mix. The water from the public mains can affect the correct setting process of the gypsum, making it faster or slower

2.- SAFETY CAUTIONS:

Do not inhale the powder. It's recommended the use of a protective mask.

Keep away from children.

Handle the product in a well-ventilated place.

The handling of this product has to be done only by qualified people

In the event of ingestion, consult a doctor.

For further information, see the safety data sheet for the product.

Storage and expiry

The stability of the product is 18-36 months depending on the package and respecting the storage conditions. (Note: The expiry date is considered with the product in his original packing and before been opened).

Storage: Keep the packages tightly closed in a dry place away from moisture.

After each use the package must be closed to prevent the moisture absorption.

Supply

Colour: vanilla

1142-201 Pail: 20 Kg

1142-261 Paper bag: 25 Kg

1142-M11 Test pack: 100 gr

Additional information

The instructions that appear on this sheet are based on our current knowledge and are given only as information without that implies a guarantee from Protechno.

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