LABOSIL

Descripción

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LABOSIL

Description

LABOSIL is a Lab Putty silicon which cures at room temperature by a condensation reaction

Applications

LABOSIL satisfies the technical requirements of the dental lab in the following applications:

- Isolation of denture/die in a flask
- Matrix for use in making temporary crowns/bridges
- Matrix for fix prostheses
- Matrix occlusal registration in articulator
- Checkbite for CAD / CAM scanning

Advantages

- Plastic consistency , easy to model
- High breaking resistance
- Easy demoulding with resins / die
- Low linear shrinkage
- High heat resistance
- Excellent reproduction of details

Application method

Before catalysis, LABOSIL is viscous and white color and the catalyst is a red gelatin paste

Add 100 parts of LABOSIL to 3 parts of CATALYST

Practically, dose 1 measuring spoon of LABOSIL for 4 cm of paste CATALYST

Use the provided measuring spoon to dose the required quantity of LABOSIL.

Knead the putty with your fingers before adding the catalyst to obtain a more plastic consistency

Add the catalyst in the recommended proportion and mix carefully by hand to get a uniform color (about 30 sec).

Then the product is ready to be used in any of the previous applications

Physical characteristics

Mixing time:	30 sec.
Working time:	1'30" approx. (ISO 4823)
Setting time:	8 min approx.
Demoulding time:	25 min approx.
Shore A hardness (24h)	78 - 82 LABOSIL 80 88 - 92 LABOSIL 90

The given times are with a room temperature between 20-22° C
An excess of catalyst or a higher temperature can speed up the setting time.

Low temperatures or insufficient catalyst can slow down the setting time.

LABOSIL silicone does not contain any kind of pigment, therefore it may have slight differences in colour between batches due to raw materials. This fact does not implies a change in the quality of the product

Presentation

LABOSIL with Catalyst

shore 80	shore 90	
5020-020	5021-020	2 Kg
5020-050	5021-050	5 Kg
5020-100	5021-100	10 Kg
5020-M20	5021-M20	Test pack

Catalyst

5010-400	40 g
----------	------

Storage and expiry

The product is stable for 24 months after manufacturing, following the storage conditions.

Batch number and expiry date can be found in a label on each package.

Keep the product in its original package and in a dry place .
Protect from moisture and sun light . Maintain under 30°C

Guarantee:

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

Our products are subject to continuous development. We therefore reserve the right to make modifications.

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