LABOSIL-2

Descripción

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LABOSIL -2

Description

LABOSIL-2 is a Lab Putty silicon which polymerizes at room temperature by an addition reaction

Applications

LABOSIL-2 has been specially designed to satisfy the technical requirements of the dental lab in the following applications:

- Isolation of denture/die in a flask
- Matrix for fix and partial prostheses
- Matrix occlusal registration in articulator
- Matrix for resin in combined and frame works prostheses
- Checkbite for CAD / CAM scanning

Advantages

- Easy to use and dose : 1 + 1
- CAD / CAM scannable in laser or white light scanners
- Plastic consistency and easy to model which allows a high duplicating precision
- Easy demoulding and cutting
- Low linear shrinkage
- High heat resistance
- No toxic: it's Platinum catalyzed

Application method

Take equal amounts of both components Catalyst (white) and Base (green) with the provided dosing spoons.
Respect the dosage 1:1

Do not exchange the dosing spoons nor the lids

Knead the putty with your fingers to get a uniform light blue color (about 30 sec).

Then the product is ready to be used during about 2 minutes

It can be demoulded in 6 minutes

Warnings

Close the packages tightly after use.

Do not exchange the dosing spoons nor the lids

Use **vinyl gloves** not latex because they may extend the setting time

Physical characteristics

Mixing time:	30 sec.
Working time:	2''
Setting time:	6'
Shore A hardness:	85
Density (23°C)	1.56 g/cm ³
Linear change variation	< 0.2 % (24 h)

The given times are with a room temperature of 23° C and an exact mixing ratio 1:1

A higher temperature can speed up the setting time.
Lower temperatures can slow down the reaction.
A higher dose of Catalyst does not accelerate the setting

Presentation

LABOSIL -2

5023-20T 2 Kg: 1 Kg Cat. + 1 Kg Base + spoons

5023-M20 Kit INTRO : 150 g Cat + 150 g Base

Storage and expiry

The product is stable for 20 months after manufacturing, following the storage conditions.

Batch number and expiry date can be found in a label on each package.

Keep the product in its original package and in a dry place.

Protect from humidity and sun light

Maintain the temperature under 30°C

Guarantee:

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet which can be requested from your usual supplier.

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