TECHNOJEL^â

Descripción

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TECHNOJEL^â

Description

Duplicating material Agar-Agar based, reversible and of one component specially suitable for investments.

Use

Suitable for making phosphate – investment duplicates in the technique of reproduction of models.
Do not use to duplicate plaster.

Advantages

Good characteristics of hardness, elasticity and stability.
Duplicates with high reliability and with great precision and definition of details.
Its great elasticity makes possible to achieve an extremely smooth model surface
Hardening without bubbles.
Reduction of costs for their use continuous and bigger output.
Easy to use.

Handling guidelines

Reproduction:

- Cut the gelatine in small pieces (about 4 cm) and introduce in the duplicating unit.
- Fuse the product up to 90°C without surpassing this temperature.
- Leave cool the product to room temperature up to 45°C (best temperature of use), instant in that the product is ready to be emptied in muffle.
- If the master model is made of gypsum, it's better to wet it to improve the emptied.
- Once emptied, the gelification time is approximately of 30 minutes in metallic muffle.
- In plastic muffle this time is prolonged substantially.

Duplicating:

- Once hardened, cast out the model carefully
- Wash the duplicate with water and drain it well (if there remains water or other particles, the duplicate surface can present roughness
- If the gelatine is too cold, it can provoke a longer setting time and therefore a weaker and rougher duplicate.
- Fill in with the investment and vibrate for 2 minutes.

Re-use:

This product can be re-used following this indications:

- Wash the used gelatine carefully removing all investment rests
- Cut in small pieces and fuse as it was a new product
- If the gelatine it's too thick, you can add small quantities of distillate water (lost by evaporation)

When the gelatine has been re-used several times, it can lose its characteristics. Then, eliminate it and use a brand new product.

Physical characteristics

Melting point:	90 °C.	(194°F)
Solidification point:	38 °C app.	(82°F)
Temperature of use:	45 °C	(115°F)
Gelification time:	30 min.	

Presentation

5510-060 Can: 6kg.
5510-010 Can: 1kg.

Expiration and Storage

The stability of the product is of 2 years (respecting the conservation instructions.)
The product should be conserved inside its packing properly closed and in a dry atmosphere .

Additional information

The information given in this sheet is based on our level of knowledge at the moment at the date of printing and is not intended to be a warranty for the product.

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