TECHNOGOLD

Descripción

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TECHNOGOLD

Description

Low fusion investment special for casting noble alloys in crown and bridge technique.

Instructions of use

MIXING RATIO

31-34 ml water to 100 g TECHNOGOLD

MIXING

Pour water into mixing bowl, add TECHNOGOLD
Incorporate with a hand spatula 15 seconds, then mix under vacuum for 30 seconds.

INVESTING

TECHNOGOLD flows easily under light vibration. Overfilling the ring and removal of overlap after hardening improves air permeability of surface.

Allow investment to set for a minimum of 1 hour (small cylinder) before heating.

BURNOUT

Place flask in a cold furnace, heat to 300°C for 1 hour 30 minutes, then raise to 720°C for a minimum of 1 hour.

Casting must take place within 15 minutes after .

REMOVING MOULDS

Let the casting cool down to room temperature. Never use water to cool while the metal is red hot.

Packaging

Code

6670-200

Pail : 20 Kg (44 lb)

Physical characteristics

Mixing ratio powder : water 100g : 31-34ml

Mixing time (vacuum) 30 sec.

Working time 6 - 8 min.

Setting time 30 min.

Casting temperature 720°C

Total expansion max. 1.84%

Storage and expiry

The stability of the product is 18 months

Keep the packages tightly closed in a dry place.

Safety information

Avoid the formation and breathing of the dust. (it contains silica)
Work in well-ventilated areas or with aspiration

Guarantee:

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of TECHNOGOLD that can be requested from your usual supplier.

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