# ALCOVEST ALCOVEST RAPID

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	4
ISTRUZIONI PER L'USO	IT	6
NOTICE D'EMPLOI	FR	8
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		10

Consulte nuestra web  
Visit our website

Alcovest

Alcovest Rapid

# ALCOVEST/ALCOVEST RAPID

## DESCRIPCIÓN

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# ALCOVEST/ALCOVEST RAPID

## DESCRIPTION

Special polysilicate and alcohol base investment for Chrome Cobalt casting on moulds. According to ISO 15912:2016, the classification of the product is Type 2 Class 1.

This material contains quartz. **DO NOT INHALE THE POWDER!**

## INSTRUCTIONS OF USE

### PREPARING THE LIQUID

Add the contents of the small bottle (B) to the large one (A). Shake it for 1 minute and **leave the mixture to rest for 4 hours**, after which time will be ready to use. The mixed liquids will begin to react after 10 minutes. A increase of the temperature must be observed.

Storage: If you do not use all the liquid you have prepared within 10 days, we recommend storing it in the refrigerator, otherwise it might crystallize.

**WARNING: DO NOT USE CONTAINERS THAT HAVE BEEN IN CONTACT WITH PLASTER.**

### MIX PROPORTION

#### For the duplicate model:

Use 200 g of investment powder and mix with the following liquid:

15-16 ml liquid / 100 g powder

#### For the cylinder:

Use 400 g of investment powder and mix with the following liquid:

17-18 ml liquid / 100 g powder

In both cases, if the mixture is made by hand, mix it with a vibrator. Put the liquid in first, adding the powder afterwards and beating for 60 seconds.

### WORKING TIME

#### **ALCOVEST: 8-10 minutes (20-22 °C) / ALCOVEST-RAPID: 6 minutes**

The working time can be reduced when the liquid is at a higher temperature, at temperatures over 22°C, and when using contaminated utensils. The working time can be extended when the liquid and/or powder are at temperatures below 20°C.

### DUPLICATION OF MODELS

The duplication can be carried out in gelatin (TECHNOJEL®) or silicone (DUPLIFLEX).

#### FILLING THE MODELS.

Pour the mixed investment into the duplicate under vibration. To guarantee a smooth and hard surface of the models **it's important to maintain the vibration for 10 minutes.**

#### SETTING TIME

**ALCOVEST: 35-40 minutes / ALCOVEST-RAPID: 30-35 minutes**, both for models as well as cylinders.

Once the model is dry, cut off the top layer where the superfluous liquid had collected, to prevent the model from cracking when it is heated.

#### DRYING THE MODELS.

**40 minutes** in a furnace at **220-240°C.**

#### MODEL HARDENING

After drying the models they may be hardened using any system. We recommend our hardener TECHNODUR.

## FILLING THE CYLINDERS


Before filling the cylinder, paint the wax mould with a fine investment (TECHNOPAINT), or apply a wax degreasing varnish. To fill the cylinder, prepare a slightly more liquid investment than for the duplicates, to make the filling operation easier. Fill with a little vibration and **keep it for 5-10 min.**

**Setting time: ALCOVEST: 40-45 minutes / ALCOVEST-RAPID: 30-35 minutes**

Once the cylinder is dry, cut off the top layer where the superfluous liquid had collected, to prevent the cylinder from cracking when it is heated.

## HEATING

### DRYING

After drying and cutting, place the cylinder in the cold furnace, dry it at 300°C for 30-60 minutes depending on the size of the cylinder.


### HEATING

Increase the temperature by between 7 and 9 degrees per minute until it reaches:

**950°C for partial plates**

**1050°C for complete plates**

Maintain the final temperature 30-60 minutes depending on the size and number of cylinders.

If the indicated drying times and pre-heating temperatures are not observed, slight fissures in the muffle surface may occur; they are of no importance for the cast object.

## AFTER CASTING

After setting, let the cylinder cool at room temperature.

**DO NOT IMMERSE IN COLD WATER!!**

## PRECAUTION

The proportion in the powder+liquid kit is thought to an easier use for the costumer. In the practical use in the laboratory, it's possible that you need more powder or liquid. To solve this, replacement for both liquid and powder are sold.

## STORAGE AND EXPIRY

Batch number and expiry date are marked on each package.

### POWDER ALCOVEST

Store at 18-25°C.

Keep the packages tightly closed in a dry place.

The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.

### LIQUIDS A & B

Store at 7-20°C.

Keep away from sources of ignition.

Keep the bottles tightly closed after each use.

Once A&B liquids have been mixed, if prepared liquid is not used in the following 10 days, store in the refrigerator, otherwise it might crystallise.

Highly flammable liquid, toxic if inhaled, in contact with skin and if swallowed, causes damage to organs.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of ALCOVEST/ALCOVEST-RAPID. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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