
CALIBRA 3D®

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	6
ISTRUZIONI PER L'USO	IT	10
NOTICE D'EMPLOI	FR	14
INFORMACIÓN ADICIONAL		18
ADDITIONAL INFORMATION		
INFORMAZIONI SUPPLEMENTARI		
INFORMATIONS COMPLÉMENTAIRES		

Consulte nuestra web
Visit our website

CALIBRA 3D®

DESCRIPCIÓN

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CALIBRA 3D®

DESCRIPTION

CALIBRA 3D is a universal investment for the casting of 3D printed dental structures in the FAST AND NORMAL HEATING technique.

The investment is designed for the casting of crowns and bridges as well as frameworks.

The expansion can be adjusted for each alloy type as indicated in the ratio table.

The investment complies with the requirements of ISO 15912:2016.

The investment contains quartz. **DO NOT BREATHE DUST!**

INSTRUCTIONS OF USE

As with all phosphatic investments, the temperature of storage and preparation of the investment and of the mixing liquid is a determining factor in the expansion during setting, and therefore in the precision of fit and the cast surface.

The optimum preparation temperature is 20-22°C. Any variation in the optimum working temperature, whether of the powder, liquid or mixing utensils, or significant variations in the ambient temperature, affect the behaviour of the product.

LIQUID PREPARATION:

IMPORTANT:


CALIBRA 3D® must always be used with its special mixing liquid **EXPANSOR** pure or diluted.

The dilutions can be made by adding distilled or demineralised water, shaking slightly and then applying to the bottle an identifying label of the concentration made.

DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN (BECOME CRYSTALLISED).

Concentrations of the mixing liquid:

Crowns and bridges:

TYPE OF ALLOY	EXPANSOR CONCENTRATION
Precious alloys	50%
Semi-precious alloys	55%
Palladium alloys	67%
Non-precious alloys	100%

Frameworks:

TYPE OF ALLOY	EXPANSOR CONCENTRATION
Cr-Co	50%

PREPARATION OF THE MIXTURE

- To obtain good results, it's essential to use the precise measurements.
- For greater ease of use, CALIBRA 3D® is supplied in pre-dosed sachets of 160g or 400g, and a dosing tube is supplied for the liquid.
- It's always recommended the use of the whole sachet and the measure of the liquid with the original dosing tube.
- When non-pre-dosed portions are used, weigh the exact content with a scale.

The powder/liquid proportion is 100g powder/24 ml liquid:

CYLINDER SIZE	POWDER	MIXING LIQUID
3 x	1 bag / 160 g	38.4 ml
6 x	2 bag / 160 g	76.8 ml
Frameworks	1 bag / 400 g	96 ml

WARNING: DO NOT USE BOWLS NOR CONTAINERS WHICH HAVE BEEN IN CONTACT WITH PLASTER.

MIXING

- Prepare the liquid EXPANSOR according to the concentration required (as indicated above).
- Add the powder to the liquid and mix vigorously by hand for 20-30 seconds to obtain a totally moist and uniform mass.
- Then begin with the mechanical **vacuum mixing for 60 seconds**.

FILLING THE CYLINDER

Place the printed objects with the casting sprues and spouts in a silicone cylinder. After mixing, fill the cylinders with a slight vibration, for 10-15 seconds, then stop.

Crowns and bridges: We recommend applying a thin layer of TENSIO-WAX to the 3D printed structure before investing.

Frameworks:

To avoid deformation, we recommend adding a sublingual reinforcing bar to the designs.

Preferably use fine investments for a first layer (TECHNOPAINT).

WORKING TIME

4-5 minutes from the moment that powder and liquid come in contact.

(Room temperature, powder and liquid from 20 to 22°C).

Cooler liquids prolong the setting time.




CYLINDER HEATING

PLACING IN THE FURNACE

WARNING:

Abrupt ignition of the 3D resin can lead to breakage. Place the cylinders inclined or elevated on investment plates.


- Few minutes before placing in the furnace scrape with a knife the upper surface of the cylinder, to allow the gases to escape.
- Place the cylinders in the furnace with the spouts downwards.
- Avoid contact of the cylinders with the base and walls of the oven.
- In the case of several cylinders, also avoid contact between them.
- Due to the aggressive elimination of the 3D resin, do not open the furnace while heating.

HEATING

Heating can be carried out with the normal heating (from a cold furnace) or fast heating (straight to final temperature).


FOR BIG STRUCTURES, WE RECOMMEND NORMAL HEATING.

Fast heating:


- Allow the cylinders to set for 20 minutes, scrape the top surface and place them in the preheated furnace at final casting temperature.

Normal heating:

- Once the cylinders have set at least 30 minutes, scrape the top surface and place them in the cold oven.
- Dry the cylinders at 250°C for 20-30 minutes.
- Raise to the casting temperature according to the table.

Maintain the cylinders at the desired final temperature for 30-70 minutes depending on the number and size of cylinders:

CYLINDER SIZE	HEATING TIME
3 x	30 min
6 x	50 min
9 x	70 min

Crowns and bridges:

700-750°C	For precious metals
800-850°C	For precious metal-ceramic alloys
850-900°C	For non-precious alloys

Frameworks:

950-1000°C	Cr-Co alloy
------------	-------------

DEMOULDING

After the casting has been completed, allow the cylinder to cool at room temperature. Never cool it in cold water when the metal is red.

STORAGE AND EXPIRY

Batch number and expiry date are marked on each package. The stability of the product is 36 months (24 months for the liquid), in compliance with the storage instructions.

CALIBRA 3D® POWDER

Store at 18-25°C.

Keep the packages tightly closed in a dry place.

The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.


EXPANSOR LIQUID

Store at 18-25°C.

Prolonged exposure to temperatures lower than 5°C or higher than 35°C may produce freezing and/or destabilisation of the product.

Do not use the liquid if it shows a sign of cristallization.

Keep the bottles tightly closed after each use.

SAFETY INSTRUCTIONS

See document Safety instructions included in the packaging.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of CALIBRA 3D® and EXPANSOR. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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