
CALIBRA-EXPRESS®

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	6
ISTRUZIONI PER L'USO	IT	10
NOTICE D'EMPLOI	FR	14
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		18

Consulte nuestra web
Visit our website

CALIBRA-EXPRESS®

DESCRIPCIÓN

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CALIBRA-EXPRESS®

DESCRIPTION

CALIBRA-EXPRESS® is an ultra fine phosphate bound investment for the casting of all types of crown and bridges alloys for the fast and normal heating technique.

According to ISO 15912:2016, the classification of the product is Type 1 Class 2.

It can be also used in all-ceramics crowns.

This material contains quartz. **DO NOT INHALE THE POWDER!**

INSTRUCTIONS OF USE

As with all phosphatic investments, the temperature of storage and preparation of the investment and of the mixing liquid is a determining factor in the expansion during setting, and therefore in the precision of fit and the cast surface.

The optimum preparation temperature is 20-22°C. Any variation in the optimum working temperature, whether of the powder, liquid or mixing utensils, or significant variations in the ambient temperature, affect the behaviour of the product.

LIQUID PREPARATION:


IMPORTANT:

must always be used with its special mixing liquid **EXPANSOR**, pure or diluted.

The dilutions can be effected by adding distilled or demineralised water, shaking slightly and then applying to the bottle an identifying label of the concentration made.

DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN (BECOME CRYSTALLISED).

Concentrations of the mixing liquid:

TYPE OF ALLOY		EXPANSOR CONCENTRATION
Precious		50-70 %
Semi precious		70-80 %
Non Semi precious	Ni-Cr	80-100 %
	Cr-Co	80-100 %

PREPARATION OF THE CYLINDER

Metal ring and plastic or silicone cylinders can be used.

For 9x cylinders, use always metal ring.

- WITH metal ring:
 - Cover the inner part of the cylinder with a quartz band (for example RING-O-LINER).
(1 layer for 1x and 3x cylinders / 2 layers for 6x and 9x)
 - Moisten slightly with water the band to prevent the absorption of the investment.
- WITHOUT ring (plastic or silicone cylinder, for example RING-O-EXPAND):
Apply a fine layer of Vaseline to facilitate subsequent extraction.

MIXING PROPORTION

- To obtain good results, it's essential to use the precise measurements.
- For greater ease of use, CALIBRA-EXPRESS® is supplied in pre-dosed sachets of 160 g, and a dosing tube is supplied for the liquid .
- It's always recommended the use of the whole sachet and the measure of the liquid with the original dosing tube.
- When non-pre-dosed portions are used, weigh the exact content with a scale.

Powder / liquid proportion: 100g powder: 22 - 23 ml liquid:

CYLINDER SIZE	POWDER	MIXING LIQUID
1 x	60 g	13-14 ml
3 x	1 bag / 160 g	35-37 ml
6 x	2 bag / 160 g	71-74 ml

WARNING: DO NOT USE BOWLS NEITHER CONTAINERS WHICH HAVE BEEN IN CONTACT WITH PLASTER.

MIXING

- Prepare the liquid EXPANSOR according to the concentration required (as indicated before).
- Add the powder to the liquid and mix vigorously by hand for 15-20 seconds to obtain a totally moist and uniform mass.
- Then begin with the mechanical **vacuum mixing for 60 seconds.**

FILLING TH CYLINDER

We recommend the application of a fine layer of TENSIO-WAX (surface tension reducer) over the wax pattern before investing.

Once the mixing time is completed, fill the cylinder, vibrating gently, for 10-15 seconds, then stop.

WORKING TIME

6-7 minutes (20-22 °C) from the moment that powder and liquid contact.

- The working time can be reduced when the liquid is at a higher temperature
- With ambient temperatures above 22 °C.
- With the use of contaminated utensils.




HEATING

PLACING IN THE FURNACE

WARNING:

Wax ignition can lead to breakings. Place the cylinders inclined or over investment plates. Never place them directly on the oven base.

- Few minutes before placing in the furnace scrape with a knife the upper surface of the cylinder, to allow the gases to escape.
- Place the cylinders in the furnace with the spouts downwards.
- In the case of several cylinders, also avoid contact between them.
- Due to the aggressive elimination of the wax, do not open the furnace while heating.


HEATING

Heating can be carried out in two different ways:


Fast heating:

- After 25 – 30 minutes from the filling of the cylinders, place them in the pre-heated oven at final temperature (900°C maximum).

Normal heating:

- Place the cylinder in the furnace after the setting.
- Dry at 250°C for 40 minutes.
- Heat up at 7°C / min. to the casting temperature.

FINAL TEMPERATURES

700-750°C	For precious alloys
800-850°C	For metal-ceramic alloys
850-900°C	For non-precious alloys

Keep the cylinder at the final temperature for:

CYLINDER SIZE	TIME AT FINAL TEMPERATURE
1 x	25 min
3 x	35 min
6 x	45 min
9 x	60 min

After the casting has been completed, bench-cool the ring to room temperature.
Do not quench in cold water!!

Recommended EXPANSOR concentrations:

- For ceramics heated over the investment

Inlays, onlays, facing and crowns:	15-20 %
------------------------------------	---------

- For injected ceramic:

(We recommend the use of RING-O-EXPAND Press cylinders)

Inlays and onlays:	50-55 %
Facing:	60 %
Individual crowns:	60-65 %
Bridges:	70 %

Recommended heating time at final temperature:

CYLINDER SIZE	TIME
1 x	50 min
3 x	60 min

When more than 1 cylinder is heated together, add 10 minutes more for each cylinder.

Divesting: Wait until the cylinder cool down to room temperature and then follow the ceramic manufacturer instructions

CALIBRA-EXPRESS® can be used with all type of ceramics.

STORAGE AND EXPIRY

Batch number and expiry date are marked on each package.

CALIBRA-EXPRESS® POWDER

Store at 18-25°C.

Keep the packages tightly closed in a dry place.

The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.


EXPANSOR LIQUID

Store at 18-25°C.

Prolonged exposure to temperatures lower than 5°C or higher than 35°C may produce freezing and/or destabilisation of the product.

Do not use the liquid if it shows a sign of cristallization.

Keep the bottles tightly closed after each use. Prolonged exposure to the air may cause crystallisation.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of CALIBRA-EXPRESS® and EXPANSOR. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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