
CALIBRA-M® RAPID

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	4
ISTRUZIONI PER L'USO	IT	6
NOTICE D'EMPLOI	FR	8
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		10

CALIBRA-M® RAPID

DESCRIPCIÓN

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CALIBRA-M® RAPID

DESCRIPTION

Extremely fine phosphate-based precision investment for casting techniques of Cobalt Chrome on model and implant combined work, for both fast and conventional heating. According to ISO 15912:2016, the classification of the product is Type 2 Class 2.

This material contains quartz. **DO NOT INHALE THE POWDER!**

CHARACTERISTICS

- Reproduction of the smallest details
- Easy moulds removing
- Easy blasting
- Fast or conventional heating

INSTRUCTIONS OF USE

This investment is designed to be used with a mix of liquid EXPANSOR and distilled water at 50 % whether making **the duplicating model or the cylinder.**

DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN (BECOME CRYSTALLISED).

MIX PROPORTION

For the duplicate model:

Use 200 g of investment powder and mix with the following liquid:

17 ml EXPANSOR + 17 ml distilled water

For the cylinder:

Use 400 g of investment powder and mix with the following liquid:

34 ml EXPANSOR + 34 ml distilled water

Put the liquid in first, adding the powder afterwards and beating manually for 20-30 seconds and then mixing under vacuum for 60-90 seconds.

WARNING: DO NOT USE CONTAINERS THAT HAVE BEEN IN CONTACT WITH PLASTER.

WORKING TIME

4 or 5 minutes.

Room temperature, powder & liquid at 20-22°C.

If the liquid or gelatines are colder the setting time is longer.

DUPLICATION OF MODELS

The duplication can be carried out in gelatin (TECHNOJEL®) or silicone (DUPLIFLEX).

FILLING THE MODELS

Pour the mixed investment into the duplicate under vibration. To guarantee a smooth and hard surface of the models **it's important to maintain the vibration for 2 minutes.**

SETTING TIME

20 minutes, both for models as well as cylinders.

DRYING THE MODELS

20 minutes in a furnace at **220-240°C.**

MODEL HARDENING


After drying the models they may be hardened using any system. We recommend our hardener TECHNODUR. If silicon moulds are used, no hardening is necessary.

- Preferably use fine investments for the first coat (TECHNOPAINT) or degrease the wax sufficiently (TENSIO-WAX).
- Remember that you have to use the liquid EXPANSOR diluted with distilled water at 50%. If this recommendation is not followed, some fissures may appear in the surface of the cylinder.
- Fill in the cylinder with a low vibration.
- The minimum thickness of the cylinder wall should be 10mm (distance between the outlying clasp and the edge of the cylinder).
- Minutes before placing in the oven, it is recommended to scrape the top surface of the cylinder with a knife to facilitate the escape of gases.


HEATING

Heating can be carried out in two different ways:


Fast heating:

- Allow the cylinder to set for 20 minutes and then place the cylinder in the oven preheated at 900°C.
- **CAUTION: DUE TO THE AGGRESSIVE WAX REMOVAL, DO NOT OPEN THE OVEN DURING THE FIRST 20 MINUTES OF HEATING.**
- Maintain the final temperature of 950-1000°C for at least 45-60 minutes.

Conventional heating:

- Place the cylinders in the oven, preferably cold, after at least 30 minutes of setting.
- Keep it at 250°C for 20-30 minutes.
- Set the casting temperature between 950-1000°C.
- Maintain the final temperature for a minimum of 45-60 minutes.
- The best expansion is achieved with a heating of 5°C per minute.

If the indicated drying times and pre-heating temperatures are not observed, slight fissures in the muffle surface may occur; although this will not affect the casted piece.

AFTER CASTING

Wait until the cast metal gets to black. At this point, if desired, the cooling can be accelerated with water. Act carefully: firstly, wet slowly and when the cylinder is colder, then submerge it completely in water. This makes easier the removing of the moulds.

STORAGE AND EXPIRY

Batch number and expiry date are marked on each package.


POWDER CALIBRA-M® RAPID

Store at 18-25°C.

Keep the packages tightly closed in a dry place.

The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.

EXPANSOR LIQUID

Store at 18-25°C.

Prolonged exposure to temperatures lower than 5°C or higher than 35°C may produce freezing and/or destabilisation of the product.

Do not use the liquid if it shows a sign of cristallization.

Keep the bottles tightly closed after each use.

Prolonged exposure to the air may cause crystallisation.

SAFETY INSTRUCTIONS

See document Safety instructions included in the packaging.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of CALIBRA-M® RAPID and EXPANSOR. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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