
CALIBRA-M®

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	4
ISTRUZIONI PER L'USO	IT	6
NOTICE D'EMPLOI	FR	8
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		10

Consulte nuestra web
Visit our website

CALIBRA-M®

DESCRIPCIÓN

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CALIBRA-M®

DESCRIPTION

Extremely fine phosphate-based precision investment for casting techniques of Cobalt Chrome on model and implant combined work. According to ISO 15912:2016, the classification of the product is Type 2 Class 1.

This material contains quartz. **DO NOT INHALE THE POWDER!**

CHARACTERISTICS

- Reproduction of the smallest details
- Easy moulds removing
- Easy blasting

INSTRUCTIONS OF USE

This investment is designed to be used with a mix of liquid EXPANSOR and distilled water at 50 % whether making **the duplicating model or the cylinder.**

DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN (BECOME CRYSTALLISED).

MIX PROPORTION

For the duplicate model:

Use 200 g of investment powder and mix with the following liquid:

14-15 ml EXPANSOR + 14-15 ml distilled water

For the cylinder:

Use 400 g of investment powder and mix with the following liquid:

30-32 ml EXPANSOR + 30-32 ml distilled water

In both cases, if the mixture is made by hand, mix it with a vibrator. Put the liquid in first, adding the powder afterwards and beating for 60 seconds. If mixing under vacuum, mix for 30 sec.

WARNING: DO NOT USE CONTAINERS THAT HAVE BEEN IN CONTACT WITH PLASTER.

WORKING TIME

4 or 5 minutes.

Room temperature, powder & liquid at 20-22°C.

If the liquid or gelatines are colder the setting time is longer.

DUPLICATION OF MODELS

The duplication can be carried out in gelatin (TECHNOJEL®) or silicone (DUPLIFLEX).

FILLING THE MODELS

Pour the mixed investment into the duplicate under vibration. To guarantee a smooth and hard surface of the models **it's important to maintain the vibration for 2 minutes.**

SETTING TIME

30-35 minutes, both for models as well as cylinders.

DRYING THE MODELS

40 minutes in a furnace at **220-240°C.**

MODEL HARDENING

After drying the models they may be hardened using any system. We recommend our hardener TECHNODUR. If silicon moulds are used, no hardening is necessary.


PREPARATION FOR INVESTING

- Preferably use fine investments for the first coat (TECHNOPAINT) or degrease the wax sufficiently (TENSIO-WAX).
- Remember that you have to use the liquid EXPANSOR diluted with distilled water at 50%.
If this recommendation is not obeyed, some fissures may appear in the surface of the cylinder.
- Fill in the cylinder with a low vibration and **demould at 35 minutes**.

HEATING


DRYING

After setting for 35 minutes put the cylinders in the furnace, preferably cold, and heat to 240 °C for 40 minutes.

HEATING

Once the cylinders are dry, programme the casting temperature. We recommend no more than 950°C. Maintain this temperature for at least 35 minutes.

The best expansion rates are obtained at a heating of 5°C per minute.

If the indicated drying times and pre-heating temperatures are not observed, slight fissures in the muffle surface may occur; they are of no importance for the cast object.

AFTER CASTING

Wait until the cast metal gets to black.

Cool the cylinders with water. Act carefully: firstly, wet slowly and when the cylinder is colder, then submerge it completely in water. This makes easier the removing of the moulds.

STORAGE AND EXPIRY


Batch number and expiry date are marked on each package.

POWDER CALIBRA-M®

Store at 18-25°C.

Keep the packages tightly closed in a dry place.

The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.

EXPANSOR

Store at 18-25°C.

Prolonged exposure to temperatures lower than 5°C or higher than 35°C may produce freezing and/or destabilisation of the product.

Do not use the liquid if it shows a sign of crystallization.

Keep the bottles tightly closed after each use.

Prolonged exposure to the air may cause crystallisation.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of CALIBRA-M® and EXPANSOR. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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