
# CALIBRA-PRESS

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	6
ISTRUZIONI PER L'USO	IT	10
NOTICE D'EMPLOI	FR	14
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		18

Consulte nuestra web  
Visit our website

# CALIBRA-PRESS

## DESCRIPCIÓN

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# CALIBRA-PRESS

## DESCRIPTION

CALIBRA-PRESS is an ultra fine, graphite-free, phosphate bound investment for casting press ceramics with the technique of NORMAL OR QUICK HEATING.

According to ISO 15912:2016, the classification of the product is Type 1 Class 2.

This material contains quartz. **DO NOT INHALE THE POWDER!**

## INSTRUCTIONS OF USE

As with all phosphatic coatings, the temperature of storage and preparation of the coating and of the mixing liquid is a determining factor in the expansion during setting, and therefore in the precision of fit and the cast surface.

The optimum preparation temperature is 20-22°C. Any variation in the optimum working temperature, whether of the powder, liquid or mixing utensils, or significant variations in the ambient temperature, affect the behavior of the product.

### LIQUID PREPARATION:


#### IMPORTANT:

CALIBRA-PRESS must be always used with the special liquid **EXPANSOR-B** and distilled water.

The dilutions can be effected by adding distilled or demineralised water, shaking slightly and then applying to the bottle an identifying label of the concentration made.

**DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN (BECOME CRYSTALLISED).**

Concentrations of the mixing liquid:

**70 % EXPANSOR-B + 30 % DISTILLED WATER**

## PREPARATION OF THE CYLINDER

Metal ring and plastic or silicone cylinders can be used.

For 9x cylinders, use always metal ring.

- WITH metal ring:
  - Cover the inner part of the cylinder with a quartz band (for example RING-O-LINER). (1 layer for 1x and 3x cylinders / 2 layers for 6x and 9x)
  - Moisten slightly with water the band to prevent the absorption of the investment.
- WITHOUT ring (plastic or silicone cylinder, for example RING-O-EXPAND Press):
  - Apply a fine layer of Vaseline to facilitate subsequent extraction.

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- To obtain good results, it is essential to use the precise measurements.
- For greater ease of use, CALIBRA-PRESS is supplied in pre-dosed sachets of 100 g, and a dosing tube is supplied for the liquid.
- It is recommended always to use whole sachets and to measure the liquid with the original dosing tube.
- In the case of using non-pre-dosed portions of powder, weigh the exact contents with a scale.

The powder-liquid proportion is: 100 g powder/25 ml liquid:

CYLINDER SIZE	POWDER	MIXING LIQUID
1 x	100 g	25 ml
3 x	200 g	50 ml
6 x	400 g	100 ml

**WARNING: DO NOT USE BOWLS NEITHER CONTAINERS WHICH HAVE BEEN IN CONTACT WITH PLASTER.**

## MIXING

- Prepare the mixing liquid (EXPANSOR-B + distilled water)
- Add the powder to the liquid and mix briskly by hand for 15-20 seconds to obtain a totally moist and uniform mass.
- Then pass to mechanical **vacuum mixing for 60 seconds** depending on the mixer.

## FILLING THE CYLINDER

CALIBRA-PRESS flows over wax surfaces without the need to paint them with degreasing agents. We recommend the application of a thin layer of TENSIO-WAX over the wax cast before the coating.

When the mixing time is completed, fill the cylinder, vibrating gently, for 10-15 seconds, then stop.

## WORKING TIME

**7-8 minutos** (20-22 °C) from the moment that powder and liquid contact

- The working time can be reduced when the liquid is at a higher temperature
- With ambient temperatures above 22 °C.
- With the use of contaminated utensils.




**HEATING**

PLACING IN THE FURNACE

**WARNING:**

Wax ignition can lead to breakings. Place the cylinders inclined or over investment plates. Never place them directly on the oven base.

- Few minutes before placing in the furnace scrape with a knife the upper surface of the cylinder, to allow the gases to escape.
- Place the cylinders in the furnace with the spouts downwards.
- In the case of several cylinders, also avoid contact between them.
- Due to the aggressive elimination of the wax, do not open the furnace while heating.

HEATING

Heating can be carried out in two different ways:

**Fast heating:**

- Placing the cylinder in the furnace between the 20 and 30 minutes after the mixing at the final casting temperature.


**Normal heating:**

- Place the cylinder in the furnace after the setting.
- Dry at 250°C for 40 minutes.
- Heat up at 7°C / min. to the casting temperature.

FINAL TEMPERATURES

850°C The casting temperature

Maintain the cylinder at the final temperature for:

CYLINDER SIZE	HEATING TIME
1 x	25 min
3 x	35 min

After the casting has been completed, bench-cool the ring to room temperature. **Do not quench in cold water!!**

**CALIBRA-PRESS POWDER**

Store at 18–25°C.

Keep the packages tightly closed in a dry place.

The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.


**EXPANSOR-B LIQUID**

Store at 18–25°C.

Prolonged exposure to temperatures lower than 5°C or higher than 35°C may produce freezing and/or destabilisation of the product.

Do not use the liquid if it shows a sign of cristallization.

Keep the bottles tightly closed after each use. Prolonged exposure to the air may cause crystallisation.

Batch number and expiry date are marked on each package.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of CALIBRA-PRESS and EXPANSOR-B. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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