
CALIBRA-PRESS

HOJA DE INSTRUCCIONES

ES

2

WORKING INSTRUCTIONS

EN

6

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IT

10

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FR

14

INFORMACIÓN ADICIONAL

18

ADDITIONAL INFORMATION

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INFORMATIONS COMPLÉMENTAIRES

CALIBRA-PRESS

DESCRIPCIÓN

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CALIBRA-PRESS

DESCRIPTION

CALIBRA-PRESS is an ultra fine, graphite-free, phosphate bound investment for casting press ceramics with the technique of NORMAL OR QUICK HEATING.

According to ISO 15912:2016, the classification of the product is Type 1 Class 2.

This material contains quartz. **DO NOT INHALE THE POWDER!**

INSTRUCTIONS OF USE

As with all phosphatic coatings, the temperature of storage and preparation of the coating and of the mixing liquid is a determining factor in the expansion during setting, and therefore in the precision of fit and the cast surface.

The optimum preparation temperature is 20-22°C. Any variation in the optimum working temperature, whether of the powder, liquid or mixing utensils, or significant variations in the ambient temperature, affect the behavior of the product.

LIQUID PREPARATION:

IMPORTANT:

CALIBRA-PRESS must be always used with the special liquid **EXPANSOR-B** and distilled water.

The dilutions can be effected by adding distilled or demineralised water, shaking slightly and then applying to the bottle an identifying label of the concentration made.

DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN (BECOME CRYSTALLISED).

Concentrations of the mixing liquid:

70 % EXPANSOR-B + 30 % DISTILLED WATER

PREPARATION OF THE CYLINDER

Metal ring and plastic or silicone cylinders can be used.

For 9x cylinders, use always metal ring.

- WITH metal ring:

- Cover the inner part of the cylinder with a quartz band (for example RING-O-LINER).

- (1 layer for 1x and 3x cylinders / 2 layers for 6x and 9x)

- Moisten slightly with water the band to prevent the absorption of the investment.

- WITHOUT ring (plastic or silicone cylinder, for example RING-O-EXPAND Press):

- Apply a fine layer of Vaseline to facilitate subsequent extraction.

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- To obtain good results, it is essential to use the precise measurements.
- For greater ease of use, CALIBRA-PRESS is supplied in pre-dosed sachets of 100 g, and a dosing tube is supplied for the liquid.
- It is recommended always to use whole sachets and to measure the liquid with the original dosing tube.
- In the case of using non-pre-dosed portions of powder, weigh the exact contents with a scale.

The powder-liquid proportion is: 100 g powder/25 ml liquid:

CYLINDER SIZE	POWDER	MIXING LIQUID
1 x	100 g	25 ml
3 x	200 g	50 ml
6 x	400 g	100 ml

WARNING: DO NOT USE BOWLS NEITHER CONTAINERS WHICH HAVE BEEN IN CONTACT WITH PLASTER.

MIXING

- Prepare the mixing liquid (EXPANSOR-B + distilled water)
- Add the powder to the liquid and mix briskly by hand for 15-20 seconds to obtain a totally moist and uniform mass.
- Then pass to mechanical **vacuum mixing for 60 seconds** depending on the mixer.

FILLING TH CYLINDER

CALIBRA-PRESS flows over wax surfaces without the need to paint them with degreasing agents. We recommend the application of a thin layer of TENSIO-WAX over the wax cast before the coating.

When the mixing time is completed, fill the cylinder, vibrating gently, for 10-15 seconds, then stop.

WORKING TIME

7-8 minutos (20-22 °C) from the moment that powder and liquid contact


- The working time can be reduced when the liquid is at a higher temperature
- With ambient temperatures above 22 °C.
- With the use of contaminated utensils.

PLACING IN THE FURNACE

**WARNING:**

Wax ignition can lead to breakings. Place the cylinders inclined or over investment plates. Never place them directly on the oven base.


- Few minutes before placing in the furnace scrape with a knife the upper surface of the cylinder, to allow the gases to escape.
- Place the cylinders in the furnace with the spouts downwards.
- In the case of several cylinders, also avoid contact between them.
- Due to the aggressive elimination of the wax, do not open the furnace while heating.

HEATING.

Heating can be carried out in two different ways:

**Fast heating:**


- Placing the cylinder in the furnace between the 20 and 30 minutes after the mixing at the final casting temperature.

**Normal heating:**

- Place the cylinder in the furnace after the setting.
- Dry at 250°C for 40 minutes.
- Heat up at 7°C / min. to the casting temperature.

FINAL TEMPERATURES

850°C The casting temperature

Maintain the cylinder at the final temperature for:

CYLINDER SIZE	HEATING TIME
1 x	25 min
3 x	35 min

After the casting has been completed, bench-cool the ring to room temperature. **Do not quench in cold water!!**

**CALIBRA-PRESS POWDER**

Store at 18-25°C.

Keep the packages tightly closed in a dry place.


The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.

EXPANSOR-B LIQUID

Store at 18-25°C.

Prolonged exposure to temperatures lower than 5°C or higher than 35°C may produce freezing and/or destabilisation of the product.

Do not use the liquid if it shows a sign of crystallization.

Keep the bottles tightly closed after each use. Prolonged exposure to the air may cause crystallisation.

Batch number and expiry date are marked on each package.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of CALIBRA-PRESS and EXPANSOR-B. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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