
CERAMVEST® HI-SPEED

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	6
ISTRUZIONI PER L'USO	IT	10
NOTICE D'EMPLOI	FR	14
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		18

Consulte nuestra web
Visit our website

CERAMVEST® HI-SPEED

DESCRIPCIÓN

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CERAMVEST® HI-SPEED

DESCRIPTION

CERAMVEST® HI-SPEED is an ultra fine phosphate bound investment for the casting of all types of crown and bridges alloys for **THE FAST AND NORMAL HEATING TECHNIQUE**.

According to ISO 15912:2016, the classification of the product is Type 1 Class 2.

This material contains quartz. **DO NOT INHALE THE POWDER!**

WORKING CONDITIONS


As with all phosphatic investments, the temperature of storage and preparation of the investment and of the mixing liquid is a determining factor in the expansion during setting, and therefore in the precision of fit and the cast surface.

The optimum preparation temperature is 20-22°C. Any variation in the optimum working temperature, whether of the powder, liquid or mixing utensils, or significant variations in the ambient temperature, affect the behavior of the product.

LIQUID PREPARATION

The expansion can be regulated by diluting the mixing liquid **EXPANSOR-B** in different concentrations in order to obtain high precision with any type of alloy used.

Concentrations of mixing liquid:

TYPE OF ALLOY		CONCENTRATION EXPANSOR - B
Precious		50-70 %
Semi precious		70-80 %
Non precious	Ni-Cr	80 %
	Cr-Co	80-100 %

The maximum expansion is obtained by using the liquid EXPANSOR-B in a pure state (100%). To obtain lower expansions, add distilled water.

Note:

Setting under pressure chamber reduces the setting expansion by approximately 10% (at 5 bar).

IMPORTANT:

CERAMVEST® HI-SPEED must always be used with its special mixing liquid EXPANSOR-B, pure or diluted.

The dilutions can be effected by adding distilled or demineralised water, shaking slightly and then applying to the bottle an identifying label of the concentration made.

DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN (BECOME CRYSTALLISED).

PREPARATION OF THE CYLINDER

With the fast heating technique, use metal ring.

- WITH metal ring:
 - Cover the inner part of the cylinder with a quartz band (for example RING-O-LINER). (1 layer for 1x & 3x cylinders / 2 layers for 6x or bigger cylinders.)
 - Moisten slightly with water the band to prevent absorption of the investment.
- WITHOUT ring (plastic or silicone cylinder, for example RING-O-EXPAND):
Apply a fine layer of vaseline to facilitate subsequent extraction.

MIXING PROPORTIONS

- To obtain good results, it is essential to use the precise measurements.
- For greater ease of use, CERAMVEST-HI SPEED is supplied in predosed sachets of 150 g, and a dosing tube is supplied for the liquid.
- It is recommended always to use whole sachets and to measure the liquid with the original dosing tube.
- In the case of using non-pre-dosed portions of powder, weigh the exact contents with a scale.

The powder liquid ratio is: 100g powder: 26 ml liquid

CYLINDER SIZE	POWDER	MIXING LIQUID
1 x	60 g	16 ml
3 x	1 bag / 150 g	39 ml
6 x	2 bags / 150 g	78 ml
9 x	3 bags / 150 g	117 ml

WARNING: DO NOT USE BOWLS NEITHER CONTAINERS WHICH HAVE BEEN IN CONTACT WITH PLASTER.

MIXING

- Prepare the liquid EXPANSOR-B according to the concentration required (as indicated above).
- Add the powder to the liquid and mix vigorously by hand for 15-20 seconds to obtain a totally moist and uniform mass.
- Then pass to mechanical vacuum mixing for 60 seconds depending on the mixer.

FILLING THE CYLINDER

- CERAMVEST® HI-SPEED flows over wax surfaces without the need to paint them with degreasing agents. We recommend the application of a thin layer of TENSIO-WAX over the wax cast before the coating.
- When the mixing time is completed, fill the cylinder, vibrating gently, for 10-15 seconds, then stop.

WORKING TIME

6-7 minutes (20-22 °C) from the moment that powder and liquid contact.
The working time can be reduced:

- When the liquid is at a higher temperature
- With ambient temperatures above 22 °C
- With the use of contaminated utensils




HEATING

PLACING IN THE FURNACE

WARNING:

Wax ignition can lead to breakings. Place the cylinders inclined or over investment plates. Never place them directly on the oven base.

- Few minutes before placing in the furnace scrape with a knife the upper surface of the cylinder, to allow the gases to escape.
- Place the cylinders in the furnace with the spouts downwards.
- In the case of several cylinders, also avoid contact between them.
- Due to the aggressive elimination of the wax, do not open the furnace for the first 15 minutes of heating.

HEATING

Heating can be carried out in two different ways:

Fast heating:



FOR FAST HEATING, USE ALWAYS METALLIC RING

- After 20 minutes from the filling of the cylinders, place them in the oven at temperature of 850°C maximum.

Normal heating:

- Place the cylinder in the furnace after the setting.
- Dry at 250°C for 40 minutes
- Heat up at 7°C / min. to the casting temperature.

FINAL TEMPERATURES

700-750°C	For precious alloys
800-850°C	For semi-precious alloys
850-900°C	For non precious alloys

Keep the cylinder at the final temperature for:

CYLINDER SIZE	TIME AT FINAL TEMPERATURE
1 x	25 min
3 x	35 min
6 x	45 min
9 x	60 min

After the casting has been completed, bench-cool the ring to room temperature.
Do not quench in cold water!!

Batch number and expiry date are marked on each package.

CERAMVEST®HI-SPEED POWDER

Store at 18-25°C.

Keep the packages tightly closed in a dry place.

The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.


EXPANSOR-B LIQUID

Store at 18-25°C.

Prolonged exposure to temperatures lower than 5°C or higher than 35°C may produce freezing and/or destabilisation of the product.

Do not use the liquid if it shows a sign of cristallization.

Keep the bottles tightly closed after each use. Prolonged exposure to the air may cause crystallisation.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of CERAMVEST® HI-SPEED and EXPANSOR-B. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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