## TECHNOVEST SOLDER

### Descripción

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## TECHNOVEST SOLDER

### Description

Special investment for high and low melting solders

### Instructions for use

Prepare the measures of powder and water according to the mix ratio.

To get the best results it's recommendable to use distilled water. Add the powder to the liquid in a clean rubber cup and spatule vigorously for 40 seconds.

Room temperatures lower than 20° C or the use of tap water will slow the setting

### PLACING

Firstly, join the pieces to solder with wax.

Over a plastic base, bed both pieces with the investment dough, getting a single block and leaving exposed the parts that have to be solded.

Once the dough has hardened ( 4 - 7 min ) cut the block to get a quick pre-dry

### SOLDERING

Fuse the wax with a blow lamp and also heat the pieces homogeneously, so that the soldering can flow through easily

If it's necessary, brush solder joints with soldering flux, following the instructions of the manufacturer.

### Physical characteristics

Mixing proportion:	13 ml water / 100 g powder
Mixing time:	40 sec ( by hand )
Working time:	2 – 3 min
Setting expansion:	0.15 %
Thermal expansion:	1.0 % ( 800 ° C )

Working temperature, powder and liquid 21° C

### Packaging

6675-050          5 Kg ( 2 x 2,5 Kg Al bag )

### Storage and expiry

The stability of the product is 36 months  
Keep the packages tightly closed in a dry place.

### Safety information

Avoid the formation and breathing of the dust. (It contains silica ) Work in well-ventilated areas or with aspiration

#### Guarantee:

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of TECHNOVEST SOLDER that can be requested from your usual supplier.

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