# ECO-CRYL COLD

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	4
ISTRUZIONI PER L'USO	IT	6
NOTICE D'EMPLOI	FR	8
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		10

CE  
1639

# ECO-CRYL COLD

## DESCRIPCIÓN

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# ECO-CRYL COLD

## DESCRIPTION

A self-curing acrylic resin designed for rapid repairing and relining of all denture plates made from acrylic material. It is also suitable for hydro-flasking partial and full dentures.

Type 2, Class 1, ISO 20795-1:2013.

## PROPERTIES

Quick curing, excellent joint with all kind of heat-curing and self-curing resins.  
Cadmium free.

## COMPOSITION

Powder / liquid cold cure resin.

Powder: 98% methyl methacrylate-copolymer

Liquid: 94% methyl methacrylate-monomer

## INSTRUCTIONS OF USE

### WORKING TIME

#### Dough Time (for a pourable paste)

45-60 sec. At 20-21°C and 21 g powder/ 18 ml liquid mix ratio.

#### Flow Time (pourable time)

2-3 minutes at 20-21°C.

#### Setting Time

15-20 minutes at 20-21°C room temperature.

### APPLICATION METHOD

#### MOULD PREPARATION

To assure a good insulate of the resin, put the gypsum mould in hot water and, once dry, paint the surface with a spacing varnish as TECHNO SIL. Set up the repair on a plaster model, and clean up the edges of the fracture.

#### SAMPLE PREPARATION

Measure out the powder and liquid in the following ratio:

**18 ml liquid : 21 g powder**

Pour the liquid in a mixing receptacle, add the powder and mix it without producing air bubbles.

Free measuring: Put the liquid in the receptacle and add the powder slowly until all the liquid absorbs the powder. Spatule for 20 seconds without producing air bubbles. The mixture has fluid and pourable consistency when it has to be poured in the model.

## APPLICATION METHOD

### APPLICATION

Transfer the mixture of the repair using a thin bladed spatula building to the thickness required. After 2-3 minutes it will reach the plastic phase, when the mould would not be handled more.

ECO-CRYL COLD can be used also under pressure; this polymerisation allows free-pore repairing.

Once hardened (20 min) it can be polished.

## FINISHING AND POLISHING

ECO-CRYL COLD can be finished and polished using usual methods. We recommend the use of the polishing pastes POLI-R and POLI-B.

## TECHNICAL INFORMATION

Residual monomer content for ECO-CRYL COLD when cured with standard method is < 3.8 % ( test 8.8 ISO 20795-1:2013).

## STORAGE AND EXPIRY

ECO-CRYL COLD powder and liquid retain all those qualities over a period of at least 2 years.

- Batch number and expiry date are marked on each package
- Store in a dry place at a temperature not exceeding 25 °C
- Keep out of sunlight

## SAFETY INSTRUCTIONS

ECO-CRYL COLD liquid contains flammable solvents and should be handled in well-ventilated non-smoking areas, away from heat and naked flames.

DO NOT INHALE THE VAPOURS.

Contact with the skin and eyes should be avoided and suitable protective clothing should be worn, i.e. overalls, gloves and goggles. In the event of contact wash immediately with copious amounts of water. A spillage of liquid should be cleaned up immediately with sand or earth.

In the event of fire use sand, earth or a suitable extinguisher.

Store in a cool place from direct sunlight, naked flames and sparks.

Flashpoint: 10° C.

### **Contraindication and side effects**

Do not use if the patient is known to be allergic to any of the ingredients in ECO-CRYL COLD. In individual cases, local allergic reactions to methacrylate materials have been reported.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet for products ECO-CRYL COLD powder and liquid. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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