ECO-CRYL PRO

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	4
ISTRUZIONI PER L'USO	IT	6
NOTICE D'EMPLOI	FR	8
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		10

CE
1639

ECO-CRYL PRO

DESCRIPCIÓN

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ECO-CRYL PRO

DESCRIPTION

Self-curing pouring acrylic resin for the production of total and partial denture.

Type 2, Class 2, ISO 20795-1:2013.

PROPERTIES

Specially designed for the manufacture of dentures with the pouring technique with gelatine or silicone. We recommend our MUFLA DE VERTIDO and the gelatine TECHNOJEL-PRO.

Cadmium free, easy handling, fluid consistency, long working time, high mechanical values, permanent colour.

Catalyst is free of tertiary amines.

Can be used with any kind of acrylic or porcelain teeth.

INSTRUCTIONS OF USE

TECHNICAL SPECIFICATIONS

Mixing proportion:

For 1 prosthesis: 25 g powder / 19 ml liquid

Use the dosers supplied with the product.

Mixing Time:	30 sec
Pourable time:	3-4 min
Curing under pressure:	30 min (45° C / 2-3 bar)

Times depend on the room temperature. A higher temperature will reduce the working time, a lower temperature will increase the pouring time.

These values are given with a room temperature of 20-22°C.

APPLICATION METHOD

MOULD PREPARATION

Prepare the wax model on the plaster mould. Fix it to the baseplate with wax or the magnet. Close the flask and fix it with the yellow rubber band. Place also the three-plug piece.

Fill in the MUFLA DE VERTIDO through the central hole with the TECHNOJEL-PRO and let it cool down about 1 hour.

Remove the model and clean the wax from the teeth with hot water and soap. Use the cleaning box supplied (is the base of the flask) to place the dental pieces in order during the wax-up. Fix the cleaning box with the short rubber band.

Once the teeth are clean, apply the needed retentions and set them back in their place.

Pierce two pouring channels with the cutting tubes

Soak the model in warm water until no air bubbles come out. Dry the surface and paint it with a separating varnish such as TECHNOSIL, taking care not to paint the teeth. Place the model back inside the flask and close with the rubber band.

APPLICATION METHOD**PREPARING THE RESIN**

Measure out the powder and liquid in the recommended ratio
Mix with a spatula for 30 seconds, and ensure that air is not trapped in the mixture.

ECO-CRYL PRO can be poured for 3-4 minutes

Pour the resin through one of the pouring channel until the resin overflows through the other.

POLYMERISATION

Place the flask in a pressure polymerizer with the holes upwards. Fill in with warm water (45°C) up to the edge of the flask .

Cure for 30 minutes with 2-3 bar of pressure

FINISHING AND POLISHING

Once the polymerization has finished, deflask the model and remove the spill-overs.

ECO-CRYL PRO can be finished and polished using usual methods.

We recommend the use of the polishing pastes POLI-R and POLI-B.

TECHNICAL INFORMATION

Residual monomer content for ECO-CRYL PRO when cured with standard method is < 4.0 % (test 8.8 ISO 20795-1:2013).

STORAGE AND EXPIRY

ECO-CRYL PRO powder and liquid retain all those qualities over a period of at least 2 years.

- Batch number and expiry date are marked on each package

- Store in a dry place at a temperature not exceeding 25 °C

- Keep out of sunlight

SAFETY INSTRUCTIONS

ECO-CRYL PRO liquid contains flammable solvents and should be handled in well-ventilated non-smoking areas, away from heat and naked flames.

DO NOT INHALE THE VAPOURS.

Contact with the skin and eyes should be avoided and suitable protective clothing should be worn, i.e. overalls, gloves and goggles. In the event of contact wash immediately with copious amounts of water. A spillage of liquid should be cleaned up immediately with sand or earth.

In the event of fire use sand, earth or a suitable extinguisher.

Store in a cool place from direct sunlight, naked flames and sparks.

Flashpoint: 10° C.

**Contraindication and side effects**

Do not use if the patient is known to be allergic to any of the ingredients in ECO-CRYL PRO.

In individual cases, local allergic reactions to methacrylate materials have been reported.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet for products ECO-CRYL PRO powder and liquid. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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