# TECHNO SIN

## Tixo

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	4
ISTRUZIONI PER L'USO	IT	6
NOTICE D'EMPLOI	FR	8
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		10

# TECHNO SIN

## Tixo

### DESCRIPCIÓN

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## Tixo

### DESCRIPTION

Self-curing acrylic resin for removable orthodontic plates.  
Suitable for the “salt and pepper” technique with a specially thixotropic performance.  
Type 1 according to ISO 20795-2:2013.

### PROPERTIES

Easy handling, chromatic stability, high mechanic values and quick process.

### COMPOSITION

Acrylic resin in liquid and powder form.

Powder:	100 % Methylmethacrylate-copolymer
---------	------------------------------------

Liquid:	> 90 % Methylmethacrylate-monomer
---------	-----------------------------------

### INSTRUCTIONS OF USE

#### APPLICATION METHOD

Soak the plaster model in warm water for about 30 minutes shake the exceeding water and dry the surface with a paper handkerchief.

Apply a layer of separating varnish as TECHNOSIL on the model and let it dry for 10 minutes.

The usual use method is “salt and pepper” technique.

Dispense the amount of liquid to be used in the supplied plastic pot and dose with the dropper. Do not introduce the dropper in the bottle to avoid contamination. Do not return the unused liquid to the bottle. A liquid contamination will lead to polymerization inside the bottle.

**APPLICATION METHOD**

Pour on the model 10-20 drops of liquid and turn it to achieve the dispersion of the liquid over the entire surface, creating a thin film. A brush may be used as well.

Sprinkle the powder uniformly over the wet model. Shake the excess of powder.

Apply more layers following the same procedure to get the desired thick.

Smooth the surface with a monomer layer.

Finish the curing in water at 40-45 °C and 2 bar of pressure for 8-10 minutes.

To obtain pore-free plates, introduce the work into the polymeriser before 6 minutes. This time is counted from the first moment that powder and liquid come into contact.

To get a better chemical union within the layers, we advise to start the moulding in a side of the model and put the screws or the retainers between the resin layers and coat again with liquid all the plate

**FINISHING AND POLISHING**

TECHNO SIN Tixo can be finished and polished using usual methods. We recommend the use of the polishing pastes POLI-R and POLI-B.

Repairs can be carried out with the same material.

**TECHNICAL INFORMATION**

Residual monomer content for TECHNO SIN Tixo when cured with standard method is < 5 % (test 8.5 ISO 20795-2:2013).

**STORAGE AND EXPIRY**

TECHNO SIN Tixo powder and liquid retain all those qualities over a period of at least 2 years.

- Batch number and expiry date are marked on each package
- Store in a dry place at a temperature not exceeding 25 °C
- Keep out of sunlight

**SAFETY INSTRUCTIONS**

TECHNO SIN Tixo liquid contains flammable solvents and should be handled in well-ventilated non-smoking areas, away from heat and naked flames.

DO NOT INHALE THE VAPOURS.

Contact with the skin and eyes should be avoided and suitable protective clothing should be worn, i.e. overalls, gloves and goggles. In the even of contact wash immediately with copious amounts of water. A spillage of liquid should be cleaned up immediately with sand or earth.

In the even of fire use sand, earth or a suitable extinguisher.

Store in a cool place from direct sunlight, naked flames and sparks.

Flashpoint: 10° C.

**Contraindication and side effects**

Do not use if the patient is known to be allergic to any of the ingredients in TECHNO SIN Tixo.

In individual cases, local allergic reactions to methacrylate materials have been reported.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet for products TECHNO SIN Tixo powder and liquid. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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