



vestige

PUTTY SOFT-FAST

Vestige Putty Soft-Fast is a hand mixing addition silicone suitable for the filling of the impression tray. It is the first "Universal impression tray material" since it has been designed for all the additive correction fluids used in the various dental impression techniques. VESTIGE PUTTY can adapt itself to any professional's methods and manual skills as well as to all the different dental impression techniques required by every specific clinical case.

Description / Use

Putty consistency vinylpolysiloxane created for the impression taking of fixed dental prosthesis, suitable for dental bridges and crowns, inlays, onlays, partial dental prosthesis, implants, position impressions, temporary C&B impressions. Balanced final hardness to allow an easy removal from the patient's oral cavity even if in presence of pronounced undercuts. Excellent elastomeric properties in terms of elastic recovery and resistance to deformation. It can be scanned. Colour: dark purple. Fast setting time.

Standard packaging REF 1S1100

VESTIGE PUTTY SOFT-FAST:
2x300ml cans (base+catalyst)

Characteristics

- Easy to mix
- Ideal consistency
- Appropriate working and setting times
- Balanced final hardness
- Dimensional stability
- Safe to disinfect
- Easy to scan

Advantages

- Pleasant to the touch, non-greasy
- Not muco-compressive
- Suitable for both single-step and two-steps impression techniques
- Easy removal from patient's oral cavity
- Accurate and unaltered in time impressions
- Safe for dentists
- Compatible with all CAD/CAM systems

TECHNICAL DATA

TECHNICAL DATA	VESTIGE PUTTY SOFT-FAST
ISO 4823 / ADA spec. 19	Type 0 Putty consistency
Mixing ratio	1:1
Mixing time (23°C)	30"
Total working time (23°C)	1'30"
Time in mouth (35°C)	2'30"
Setting time (35°C)	4'
Hardness	62±2
Linear dimensional change	< 0,20%
Elastic recovery	> 99,5%