

TABLA DE YESOS • STONES TABLE • TABLE DE PLÂTRES



DIAMANT	KIMBERLIT	RESIN-GIPS	RESIN-GIPS PLUS	MALAKIT	RUBINIT	ROCANIT 0,08	ORTOGUIX	DURGUIX	QUICK-ROCK	FLU-STONE	ORTOGUIX PLASTER	PARIS-BLANC	MUFLA-GIPS

DATOS TÉCNICOS • TECHNICAL SHEET • DONNÉES TECHNIQUES

ISO 6873 TYPE:CLASS		4	4	4	4	4	4	3	3	3	2:1	4	2:2	2:2	2:2	
Material · Material · Matériel		Sythetic	Natural	Nat+Resin	Nat+Resin	Natural	Natural	Natural	Natural	Natural	Natural	Natural	Natural	Natural	Natural	Natural
Consistencia · Consistency · Consistance (ml)		20	20	20	20	20-23	22-24	30	27-28	28-30	30	24	37-38	50	35-40	
Mezcla mecánica al vacío · Mechanical vacuum mixing · Mélange mécanique sous vide (sec)		30	30	30	30	30	30	30	30	30	-	30	30	-	-	
Mezcla manual · Manual mixing · Mélange manuel (sec)		60	60	60	-	60	60	60	60	60	30	60	60	60	60	
Tiempo de trabajo · Working time · Temps de travail (min)		8	7	7	7	7	7	5	6	5	1.5	6	7	3	3	
Fraguado final · Final setting · Final temps de prise Vicat (min)		15	13	10	10	12	12	9	13	12	3	11	13	7	8	
Expansión · Expansion · Expansion 2 h(%)		0.08	0.08	0.10	0.10	0.10	0.10	0.08	0.13	0.16	0.01	0.05	0.16	0.29	0.30	
Resistencia compresión · Compressive strength · Résistance à la compression 1 h(MPa)		60	58	50	50	53	40	32	35	32	25	40	25	10	15	

APLICACIONES · APPLICATIONS SHEET · APPLICATIONS

Prótesis fija Fixed prosthesis Prothèse fixe	CAD/CAM	✓	✓													
	Implantes · Implants · Implants	Optional	Optional	✓	✓											
	Modelo maestro · Master cast · Maître modèle	✓	✓	✓	✓	✓										
	Muñones · Dies · Dies	✓	✓	✓	✓	✓										
	Inlays-Onlays	✓	✓	✓	✓	✓										
	Puentes · Bridges · Bridges	✓	✓	✓	✓	✓	✓									
Prótesis removable Removable prosthesis Prothèse amovible	Antagonistas · Antagonists · Antagonistes (High Res.)	✓	✓			✓	✓	✓								
	Esquelético · Frameworks · Squelettique					✓	✓	Optional		Optional						
	Antagonistas · Antagonists · Antagonistes					✓	✓	✓	Optional	✓						
Ortodoncia · Orthodontics · Orthodontie	Acrílico · Acrylic · Acrylique							✓		✓						Optional
	Bases · Bases · Bases								✓				✓			
	Mufla · Flask · Moufle												Optional	✓	✓	
Otros · Other · Autres	Articulador · Articulator · Articulateur										✓					