



QUICK-ROCK

ARTICULATION / STONES



EXPANSION
0.01%



QUICK SETTING AND VERY LOW EXPANSION

- SETTING EXPANSION EXTREMELY LOW
- HARDNESS IN ACCORDANCE TO THE NEEDED REQUIREMENTS
- ADHERENT EFFECT

STONE TYPE 2 : CLASS 1 • ISO 6873

QUICK-ROCK

SPECIAL STONE FOR MOUNTING.

ADVANTAGES

- **Quick setting and very low expansion.**
- Setting expansion extremely low
- Hardness in accordance to the needed requirements

USE INSTRUCTIONS

Proporción polvo/agua:	100 g/ 30 ml
Working temperature:	20-23°C
Mechanical vacuum mixing:	20 sec.
Manual mixing:	30 sec.

PHYSICAL CHARACTERISTICS

Working time:	1,5 min.
Initial setting:	2 min.
(Tests performed according to ISO 6873:2013)	
Final setting:	3 min.
(Tests performed according to ISO 6873:2013)	
Setting expansion (2 h.):	0.01 %
Compressive strength (1 h.):	25 N / mm ²

SALES FORMAT

WHITE

1160-100	Eco-pack (5 x 2 kg)	10 kg
1160-150	Pail	15 kg
1160-K20	Test pack (1 x 2 kg)	2 g

RELATED PRODUCTS

KIMBERLIT

Extra hard high density die stone.



ISOGIPS

Separating agent plaster-plaster. EFFECTIVE AND QUICK SEPARATION.



ARTICULADOR

Semiadjustable articulator without face-bow. Made of technical plastic.



protechno

ADVANCED DENTAL PRODUCTS



Tel. +34 972 526 169
Polígono Empordà Internacional
17469 . VILAMALLA (Girona) . Spain
protechno@protechno.com
www.protechno.com