



GINGIPRO-SIL

SILICONES



MODELS WITH A NATURAL GUM SIMULATION

- EASY AND FAST USE
- HIGH STABILITY
- NATURAL GUM COLOUR



protechno

ADVANCED PRODUCTS FOR DENTAL LABS

GINGIPRO-SIL

SILICONE BASED GINGIVAL MASK FOR THE REPRODUCTION OF SOFT TISSUE CONTOURS.

APPLICATIONS

- For models with simulated gum manufacture
- For the determination of the shape of the cervical contours
- The emergence profile of the prosthesis implant

ADVANTAGES

- Addition silicone formulation in a cartridge
- Easy and fast use
- It can be extruded directly into the mould without any bubbles

GINGIPRO-SIL Rigid

Advantages

- Its optimum hardness allows trimming of the edges in the model finishing
- Easy to trim
- Accurate margin definition
- Precision and dimensional stability

GINGIPRO-SIL Elastic

Advantages

- The high elastic module facilitates model removal and reinsertion
- Excellent aesthetic result due to its translucency
- Precision and dimensional stability
- Can be used both with direct and indirect technique

SALES FORMAT

GINGIPRO-SIL Rigid	5072-K10
GINGIPRO-SIL Elastic	5072-K11
2 x 50 ml + 12 mixing tips yellow + 1 separator + 12 oral tips yellow	

Accessories

5073-000	Dispensing device 1:1 / 2:1	1 un.
5073-CA4	Mixing tips yellow	48 un.
5073-PA4	Oral tips yellow	48 un.
5072-511	Separator 10 ml	1 un.

PHYSICAL CHARACTERISTICS

GINGIPRO-SIL Rigid

Mixing ratio Base: Catalyst	1 : 1
Total working time (23°C/73°F)	1'30''
Setting time (23°C/73°F)	7-8'
Hardness (Shore A after 24 h)	70

GINGIPRO-SIL Elastic

Mixing ratio Base: Catalyst	1 : 1
Total working time (23°C/73°F)	1'30''
Setting time (23°C/73°F)	7-8'
Hardness (Shore A after 24 h)	40



ADVANCED PRODUCTS FOR DENTAL LABS



Tel. +34 972 526 169

Fax +34 972 525 586

Polígono Empordà Internacional
17469 . VILAMALLA (Girona) . Spain

protechno@protechno.com

www.protechno.com