



POLI-R

POLISHING

NEW CONTAINER

AN
ALTERNATIVE
TO PUMICE
STONE



FORGET ABOUT **POWDER** AND CHANGE TO **GEL!**

- PUMICE SUBSTITUTE
- CLEANER, FASTER, EASIER
- POLISHING AND SHINING WITH THE SAME PRODUCT



protechno

ADVANCED PRODUCTS FOR DENTAL LABS

POLI-R

**POLISHING GEL FOR ACRYLIC WORK.
PUMICE SUBSTITUTE.**

ADVANTAGES


- Fast polishing speed whitout dust
- Contains bactericide: CHLORHEXIDINE
- It is easily removed with hot water

HOW TO USE

POLI-R



[1]
ABRASIVE POLISHING
with felt cone or brush.



[2]
FINE POLISHING
with cotton brush.

Gloss finishing



POLI-B



[3]
GLOSS POLISHING
with cotton brush.

High gloss finishing



SALES FORMAT

7880-002

400 g (0.9 lb)

7880-M10

Test pack 120 g (4.2 oz)

ASSOCIATED PRODUCTS

POLI-B

Polishing paste for resins and precious alloys.



ADVANCED PRODUCTS FOR DENTAL LABS



Tel. +34 972 526 169
Polígono Empordà Internacional
17469 . VILAMALLA (Girona) . Spain
protechno@protechno.com
www.protechno.com