



protechno

by

zeiser

PERFORADORA P10
PIN-DRILL P10



PROTECHNO by ZEISER



Sistema para confeccionar modelos desmontables.

System for making demountable models.

Perforación automática de pines, de gran precisión.

High automatic precision drill machine for pin system.

Soporte para placa base electro-magnético, compatible con otras marcas.

Baseplate holder electro-magnetic, compatible with other brands.

Sistema de autofijación.

Self-fixing system.

Provista de guía láser.



Equipped with laser guide.

Rápida y de fácil manejo.

Fast and easy to use.

Colores disponibles de las placas:

Available colours of the base plates:



Soporte/Holder

CARACTERÍSTICAS TÉCNICAS

TECHNICAL DATA

Dimensiones / Dimension (cm): 25 x 18,5 x 35 h

Peso /Weight: 10 kg

Tensión / Mains supply: 230 V

Láser Clase II /Laser Class II: < 1 mW



protechno

ADVANCED PRODUCTS FOR DENTAL LABS

Tel. +34 972 526 169

Polígono Empordà Internacional
17469 . VILAMALLA (Girona) . Spain

protechno@protechno.com

www.protechno.com