

Safety data sheet
according to 91/155/EC

Date of printing: 24/10/11

Revised at: 01/12/08

1. Identification of the substance / preparation and company / undertaking

product information:

commercial name: VETAS DE FIBRA - VEINS

article Nr.: 9090

information on manufacturer / supplier

manufacturer / supplier: FAMADENT S.L.
address: Polígono Empordà Internacional
nat. registration / post code / place: E-17469 Vilamallà, Girona, Spain
telephone: +34-972 52 61 69*
fax: +34-972 52 55 86
E-mail: protechno@protechno.es

enquiry department: Technical department **tel:** +34-972 52 61 69

emergency enquiry: Technical department **tel:** +34-972 52 61 69

2. Composition properties (component substances)

CAS-Nr.	Name acc. to EC guidelines	%	TLV (units)	PEL (units)
9004-35-7	Cellulose Acetate	< 99	NE	NE

3. Hazards identification

Principal Hazard:

Primary routes of entry: Eyes or skin (no absorption) or inhalation

carcinogenicity: None of the components of this material are listed by IARC, NTP, OSHA or ACIGH as a carcinogen

effects of over exposure: Fibers are considered to be articles under the OSHA Hazard Communication Standard.

4. First Aid measures

- inhalation:** Remove to fresh air. Get medical help if discomfort persists
- ingestion:** Rinse mouth out with water. Call doctor if amount was large
- contact with the eyes:** Flush with water for 15 minutes, including under eyelids.
- contact with the skin:** Wash skin in running water and mild soap.



protechno

Poligono Empordà Internacional
17469 Vilamalla, Girona, Spain

Safety data sheet
according to 91/155/EC

Date of printing: 24/10/11

Revised at: 01/12/08

5. Fire- fighting measures

- . **extinguisher method:** Water, carbon dioxide, dry chemical
- . **unsuitable extinguisher method:** None known
- . **special fire fighting procedures:** Fire-fighters should wear self-contained breathing apparatus

6. Accidental release measures

- . **precautions with reference to persons:**
Use in well ventilated areas.
- . **precautions with reference to environment:** None
- . **procedure for cleaning / absorption:**
Use vacuum cleaner with appropriate filter or sweep up .
- . **additional information:**
Repair damaged packages carefully

7. Handling and storage

- . **handling :**
Handle carefully in a ventilated area
Keep the packages tightly closed
Wash face and hands thoroughly with soap and water after use and before eating, drinking or smoking
- . **storage:**
Store in a cool dry place.
Keep from contact with oxidizing materials

8. Exposure controls / personal protection

- . **personal protective equipment:** Measures for general protection and hygiene.
- . **respiratory protection:** Use type for particulates not otherwise classified, if needed
- . **protection of the hands:** Rubber gloves
- . **protection of the eyes:** Protective glasses (recommendable)



protechno

Poligono Empordà Internacional
17469 Vilamallà, Girona, Spain

Safety data sheet
according to 91/155/EC

Date of printing: 24/10/11

Revised at: 01/12/08

9. Physical and chemical properties

- . **form:** Veins
- . **colour:** Red
- . **odour:** Little or no odour
- . **boiling point:** not applicable
- . **thermal decomposition:** not applicable
- . **specific gravity (H₂O = 1):** 1.3
- . **solubility in water:** Non-soluble (20° C)
- . **vapour density:** not applicable
- . **pH value:** not applicable

10. Stability and reactivity

- . **stability:** The product is stable
- . **conditions to be avoided:** Heating above 240 ° C
- . **substances to be avoided:** Strong oxidizing agents
- . **decomposition products :** Oxides of carbon
- . **hazardous polymerization:** Will not occur

11. Toxicological information

Although this product may contain pigments having TLV's and/or PEL's, it is judged that none of the pigments used will be realised in form or concentrations likely to exceed exposure limits or procedure significant adverse health effects

12. Ecological information

- . **data relating to disposal:** Do not discharge into drains or the environment
- . **bio-degradability:** No data available

13. Disposal considerations

Material may be disposed of at approved landfill or incinerated under approved conditions in compliance with local regulations.

14. Transport information

Material not regulated following IMDG, IATA, ADR and UN regulations.
Store away from foodstuff



protechno

Poligono Empordà Internacional
17469 Vilamallà, Girona, Spain

Safety data sheet
according to 91/155/EC

Date of printing: 24/10/11

Revised at: 01/12/08

15. Regulatory information

No labelling is required in accordance with the 67/548/EEC directive.

- **Risk phrases:** None
- **Safety phrases:** None

16. Other information

The information is based on our level of knowledge at the moment at the date of printing and is not intended to be a warranty for the product

N.A. : Not applied
N.E. : Not established