

# GIPSLISTE • TABELLA DEI GESSI • STONES TABLE



| DIAMANT | KIMBERLIT | RESIN-GIPS | RESIN-GIPS PLUS | MALAKIT | RUBINIT | ROCANIT 0,08 | ORTOGUIX | DURGUIX  | QUICK-ROCK | FLU-STONE | ORTOGUIX PLASTER | PARIS-BLANC | MUFLA-GIPS |
|---------|-----------|------------|-----------------|---------|---------|--------------|----------|----------|------------|-----------|------------------|-------------|------------|
|         |           |            |                 |         |         |              |          | <br><br> |            | <br><br>  |                  |             |            |

## TECHNISCHE DATEN • DATI TECNICI • TECHNICAL SHEET

| ISO 6873 TYP-KLASSE   |  | 4        | 4       | 4         | 4         | 4       | 4       | 3       | 3       | 3       | 2:1     | 4       | 2:2     | 2:2     | 2:2     |         |
|---|--|----------|---------|-----------|-----------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| Material · Materiale · Material   |  | Sythetic | Natural | Nat+Resin | Nat+Resin | Natural | Natural | Natural | Natural | Natural | Natural | Natural | Natural | Natural | Natural | Natural |
| Konsistenz · Consistenza · Consistency (ml)   |  | 20       | 20      | 20        | 20        | 20-23   | 22-24   | 30      | 27-28   | 28-30   | 30      | 24      | 37-38   | 50      | 35-40   |         |
| Mechanische Vakuummischung · Miscelazione meccanica sottovuoto · Mechanical vacuum mixing (sec) |  | 30       | 30      | 30        | 30        | 30      | 30      | 30      | 30      | 30      | -       | 30      | 30      | -       | -       |         |
| Manuelle Mischung · Miscelazione manuale · Manual mixing (sec)                                  |  | 60       | 60      | 60        | -         | 60      | 60      | 60      | 60      | 60      | 30      | 60      | 60      | 60      | 60      |         |
| Verarbeitungszeit · Tempo di lavoro · Working time (min)  |  | 8        | 7       | 7         | 7         | 7       | 7       | 5       | 6       | 5       | 1.5     | 6       | 7       | 3       | 3       |         |
| Abbindungszeit · Indurimento finale · Final setting (min)                                       |  | 15       | 13      | 10        | 10        | 12      | 12      | 9       | 13      | 12      | 3       | 11      | 13      | 7       | 8       |         |
| Abbindeexpansion · Espansione · Expansion 2 h(%)  |  | 0.08     | 0.08    | 0.10      | 0.10      | 0.10    | 0.10    | 0.08    | 0.13    | 0.16    | 0.01    | 0.05    | 0.16    | 0.29    | 0.30    |         |
| Druckfestigkeit · Resistenza alla compressione · Compressive strength 1 h(MPa)                  |  | 60       | 58      | 50        | 50        | 53      | 40      | 32      | 35      | 32      | 25      | 40      | 25      | 10      | 15      |         |

## ANWENDUNGEN · APPLICAZIONI · APPLICATIONS SHEET

|  |  |          |          |          |          |          |          |          |          |          |   |   |          |   |   |          |
|--|--|----------|----------|----------|----------|----------|----------|----------|----------|----------|---|---|----------|---|---|----------|
| Fix Prothesen<br>Protesi fissa<br>Fixed prosthesis                     | CAD/CAM  | ✓        | ✓        |          |          |          |          |          |          |          |   |   |          |   |   |          |
|  | Implantate · Impianti · Implants                     | Optional | Optional | ✓        | ✓        |          |          |          |          |          |   |   |          |   |   |          |
|  | Meistermodell · Modello principale · Master cast     | ✓        | ✓        | ✓        | ✓        | ✓        |          |          |          |          |   |   |          |   |   |          |
|  | Stümpfe · Monconi · Dies                             | ✓        | ✓        | ✓        | ✓        | ✓        |          |          |          |          |   |   |          |   |   |          |
|  | Inlays-Onlays  | ✓        | ✓        | ✓        | ✓        | ✓        |          |          |          |          |   |   |          |   |   |          |
|  | Brücken · Ponti · Bridges                            | ✓        | ✓        | ✓        | ✓        | ✓        | ✓        |          |          |          |   |   |          |   |   |          |
| Herausnehmbare Prothesen<br>Protesi rimovibile<br>Removable prosthesis | Antagonisten · Antagonisti · Antagonists (Hole Res.) | ✓        | ✓        |          |          | ✓        | ✓        | ✓        |          |          |   |   |          |   |   |          |
|  | Skelett · Scheletrati · Frameworks                   |          |          |          |          | ✓        | ✓        | Optional |          | Optional |   |   |          |   |   |          |
|  | Antagonisten · Antagonisti · Antagonists             |          |          |          |          | ✓        | ✓        | ✓        | Optional | ✓        |   |   |          |   |   |          |
|  | Acryl · Acrilico · Acrylic                           |          |          |          |          |          |          | ✓        |          | ✓        |   |   |          |   |   | Optional |
| Orthodontie · Ortodonzia · Orthodontics                                |  |          |          |          |          |          | Optional |          | ✓        |          |   |   | ✓        |   |   |          |
| Andere · Altro · Other   | Basis · Base · Bases                                 |          |          | Optional | Optional | Optional | Optional | Optional |          |          |   | ✓ |          |   |   |          |
|  | Küvette · Muffola · Flask                            |          |          |          |          |          |          |          |          | Optional |   |   | Optional | ✓ | ✓ |          |
|  | Artikulator · Articolatore · Articulator             |          |          |          |          |          |          |          |          |          | ✓ |   |          |   |   |          |