



# vestige

## BITE

Vestige Bite is an impression material for bite registration. Its initial imperceptible consistency allows patients to bite naturally. VESTIGE BITE includes the characteristics the perfect material to perform correct bite registrations must have: creamy while extruding, thixotropic and stable, fast set and rigid after setting.

### Description / Use

Bite registration vinylpolysiloxane used for registering quadrants or the full-arch. It is recognizable by its rigidity, precision and setting times. Easy to use for dentists and reliable for dental technicians. Can be scanned. Colour: green. Fast Setting time.

#### Standard packaging REF 1S4000

VESTIGE BITE: 2x50 ml cartridges + 6 mixing tips

### Characteristics

- Thixotropic
- Rapid
- Rigid
- Dimensional stability
- Safe to disinfect
- Easy to scan
- Easy to trim

### Advantages

- Doesn't slip on the occlusal plane
- Avoids possible recording errors caused by micro-movements
- No thickness risks
- Accurate and unaltered in time impressions
- Safe for dentists
- Compatible with all CAD/CAM systems
- Not altered by trimming

### TECHNICAL DATA

TECHNICAL DATA	VESTIGE BITE
Mixing ratio	1:1
Mixing time (23°C)	Automix
Total working time (35°C)	30"
Time in mouth (35°C)	60"
Setting time(35°C)	1'30"
Hardness	95 Shore "A"
Linear dimensional change	< 0,05%
Details reproduction	20µm