



vestige

IMPLA CLEAR

Vestige IMPLA CLEAR is a transparent medium consistency addition silicone, suitable for direct impressions in the oral cavity and for the creation of transparent masks.

Description / Use

Vinylpolysiloxane for the making of transparent masks, occlusal key, bite rebasing, and model stabilisation in guided implantology.

Standard packaging REF 1S6100

VESTIGE IMPLA CLEAR: 2x50ml cartridges + 6 mixing tips

Characteristics

- Fluid
- Transparent
- Rigid after setting
- Dimensional stability
- Safe to disinfect
- Easy to mix

Advantages

- Flowing, it enables reproduction without incorporating bubbles
- It does not interfere with resin polymerization and reduces its inhibited layer
- It does not deform
- Accurate and unaltered in time impressions
- Safe for dentists
- Packaged in standard 50ml cartridge

TECHNICAL DATA

TECHNICAL DATA	VESTIGE IMPLA CLEAR
ISO 4823 / ADA spec. 19	Type 2 medium bodied consistency
Mixing ratio	1:1
Mixing time (23°C)	Automix
Total working time (23°C)	40"
Time in mouth (35°C)	2'
Setting time (35°C)	3'
Hardness	70 ± 2 Shore "A"
Linear dimensional change	< 0,20%
Elastic recovery	> 99,5%